#4219 EMTRAC O
Embryo Transfer Catheter

0120

Important information, read carefully before use

Disposable – discard after single use
Re-use may cause sexually transmitted infections
Caution: federal (USA) law restricts this device to sale by or on the order of a physician
Sterile: unless pouch is damaged or opened

INDICATIONS FOR USE:

Introduction of embryos into the uterine cavity.

CONTRAINDICATIONS:

Do not use in the presence of cervical infection.
Do not use in the presence or after recent Pelvic Inflammatory disease.
Do not use for intrafallopian procedures.

WARNINGS:

Not sterile if pouch is damaged or opened, discard immediately.

PRECAUTIONS:

By virtue of its flexibility and small diameter, the Emtrac O can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroflexion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

In no case should the device be forced against digitally felt resistance. If passage of the Emtrac O through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered.

The Emtrac O should be used by healthcare professionals only

ADVERSE REACTIONS:

No adverse reactions associated with the use of The #4219 Emtrac O have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

Manufactured by:
Gynetics Medical Products N.V.
Rembert Dodoensstraat 51
3920 Lommel, Belgium

HOW SUPPLIED:

#4219 Emtrac O

PACKAGED:

3 boxes of 25 individually sterile packed, disposable sets of devices.

DIRECTIONS FOR USE:

1. With the patient in the dorsal lithotomy position expose the uterine cervix to view with a vaginal speculum.
2. Having previously determined the depth, set positioner on the guiding catheter (white handle), the Emtrac O-A.
3. Having previously determined direction of the patient's uterus, put the mandrel (blue handle), the Emtrac O-B, in the guiding catheter (white handle) Emtrac O-A to pre-shape the guiding catheter.
4. Gently advance the catheter into the external cervical os and advance it through the cervical canal near the uterine cavity until resistance is felt.
5. Remove the Mandrel (blue handle) Emtrac O-B when position is reached
6. Load the embryo(s) into the loading catheter (pink handle), The #4219 Emtrac O, C catheter following your standard protocol. Insert the loading catheter (white handle). The marking on the proximal end of the tube indicates when the distal tip of the guiding catheter is reached.
7. The distal tip of the loading catheter (pink handle) should then be at the uterine fundus.

Note: If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. In no case should force be used against resistance to achieve insertion.

8. Slowly release the embryos into the uterus.
9. Following the complete expulsion of the embryos, slowly withdraw the catheter from the uterus and discard it. Remove the speculum and allow the patient to rest.
10. Once the set is removed from the patient, inspect the loading catheter (pink handle) #4219 Emtrac O C catheter to verify embryo release.

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