
**#2000 SEMTRAC W**

Soft Embryo Transfer Catheter

Important information, read carefully before use

**Disposable – discard after single use**  
**Re-use may cause sexually transmitted infections**

**Caution: federal (USA) law restricts this device to sale by or on the order of a physician.**

The trial devices are intended to support the use of an approved embryo transfer catheter. These tools are not intended to directly transfer gametes/embryos into the uterine cavity; rather, they are to aid in traversing a difficult or stenotic cervix or to perform a trial/mock transfer to better identify the uterine environment.

**Sterile: unless pouch is damaged or opened**

**INDICATIONS FOR USE:**

Introduction of embryos into the uterine cavity.

**CONTRAINDICATIONS:**

Do not use in the presence of cervical infection.

Do not use in the presence or after recent Pelvic Inflammatory disease.

Do not use for intrafallopian procedures.

**WARNINGS:**

Not sterile if pouch is damaged or opened, discard immediately.

**PRECAUTIONS:**

By virtue of its flexibility and small diameter, the #2000 Semtrac Set W can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine antelexion or retroflexion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

In no case should the device be forced against digitally felt resistance.

If passage of the #2000 Semtrac W through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered.

The Semtrac W should be used by healthcare professionals only

**ADVERSE REACTIONS:**

No adverse reactions associated with the use of the #2000 Semtrac W have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

**Manufactured by:**

Gynétics Medical Products N.V.  
 Rembert Dodoensstraat 51  
 3920 Lommel, Belgium

**HOW SUPPLIED:**

- |                              |                              |
|------------------------------|------------------------------|
| - #2000 Semtrac W-240        | - #2000 Semtrac W-180        |
| - #2000 Semtrac W-240 TIS    | - #2000 Semtrac W-180 TIS    |
| - #2000 Semtrac W-240 SP     | - #2000 Semtrac W-180 SP     |
| - #2000 Semtrac W-240 SP TIS | - #2000 Semtrac W-180 SP TIS |
| - #2000 Semtrac W-240 Trial  | - #2000 Semtrac W-180 Trial  |
| - #2000 B W-240              | - #2000 B W-180              |

**PACKAGED:**

- **#2000 Semtrac W-240 (TIS)/ #2000 Semtrac W-180 (TIS)/ #2000 Semtrac W-180 Trial/ #2000 Semtrac W-240 Trial:**  
3 boxes of 25 individually sterile packed, disposable sets of devices and 5 mandrels.
- **#2000 Semtrac W-240 SP (TIS)/ #2000 Semtrac W-180 SP (TIS):**  
3 boxes of 25 individually sterile packed, disposable sets of devices and 5 mandrels. Packaged in a double pouch.
- **#2000 B W-240/ #2000 B W-180:**  
5 mandrels per 25 #2000 Semtrac W-240/ #2000 Semtrac W-180/ #2000 Semtrac W-180 Trial/ #2000 Semtrac W-240 Trial

**DIRECTIONS FOR USE FOR THE #2000 SEMTRAC W-240 (TIS) / #2000 SEMTRAC W-240 SP (TIS) / #2000 SEMTRAC W-180 (TIS) / #2000 SEMTRAC W-180 SP (TIS):**

1. With the patient in the dorsal lithotomy position expose the uterine cervix to view with a vaginal speculum.
  2. Having previously determined the depth and direction of the patient's uterus, gently advance the guiding catheter (white handle), into the external cervical os and advance it near the cervical canal into the uterine cavity until resistance is felt.
- Note:** If necessary, the use of the mandrel (blue handle) can facilitate insertion of the catheter. This by preforming and reinforcing the guiding catheter.
3. Remove the Mandrel from the guiding catheter (if applicable)
  4. Load the embryo(s) into the loading catheter (pink handle), following your standard protocol.
  5. Insert loading catheter (pink handle) into the guiding catheter (white handle).
  6. The markings on the proximal end of the loading catheter (pink handle) indicate the distance the loading cannula has advanced beyond the tip of guiding catheter (white handle).
  7. The distal tip of the catheter should then be at the uterine fundus. **Note:** The tip of the TIS loading catheters can be viewed under ultrasound.

**Note:** If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. In no case should force be used against resistance to achieve insertion.

8. Slowly release the embryos into the uterus.
9. Following the complete expulsion of the embryos, slowly withdraw the catheter from the uterus and discard it. Remove the speculum and allow the patient to rest.
10. Once the set is removed from the patient, inspect the loading catheter (pink handle) to verify embryo release.
11. Discard after use

**DIRECTIONS FOR USE FOR THE #2000 SEMTRAC W-180 TRIAL/ #2000 SEMTRAC W-240 TRIAL:**

1. Place patient in lithotomy or knee/chest position.
2. Insert vaginal speculum to expose cervix.
3. Use malleable/flexible obturator with guiding catheter to traverse cervix if necessary. Remove malleable/flexible obturator prior to placement of trial/mock catheter. (May be performed under ultrasound guidance if desired.)
4. Introduce trial/mock catheter through guiding catheter.
5. Remove trial/mock catheter from guiding catheter.
6. If so desired, move on to embryo transfer using approved catheter and refer to instructions for use for that approved catheter.
7. Upon completion of trial/mock procedure, remove the catheter/catheters, and discard.

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