**GYNÉTICS®**  
Cu250 / Cu375

Gynetics® IUD Cu250  
Gynetics® IUD Cu375

**Information on the intrauterine contraceptive device, a method to prevent pregnancy:**

World-wide, millions of women use an intrauterine contraceptive device as a method to prevent pregnancy. There are different kinds of intrauterine contraceptive devices- or IUDs. IUD is the abbreviation for intra-uterine device. Studies show that the IUD is a reliable and safe method of anti-conception. The Gynetics® IUD is made of plastic and is partially wrapped with copper. The IUD is shaped according to the anatomy of the uterus.

**IUD picture 1: (A) side arms (B) stem with copper (C) nylon threads.**

**How does the Gynetics® IUD Cu250/ Cu375 work?**

Compared to other temporary and reversible methods, the IUD is almost as safe as the birth control pill. Although no method to prevent pregnancy is a 100% safe, a Gynetics® IUD is very reliable when it comes to protecting you against pregnancy. An IUD does not protect you from sexually transmitted diseases.

The mechanism of the IUD is still not completely clear. It is likely that the presence of the IUD disturbs the normal function of the male sperm cells, which prevents fertilisation of the egg. It is also likely that the copper effects the development of the egg, which prevents fertilisation of the egg. For this reason IUDs are not thought to induce abortion.

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**How is the Gynetics® IUD Cu250/ Cu375 inserted?**

The best time for insertion of the Gynetics® IUD is during menstruation. Most women tolerate insertion well. Before your physician inserts the IUD, the womb is examined and the length ascertained. The physician has also examined you carefully to see whether the IUD is a suitable method to prevent pregnancy for you.

The IUD is carefully inserted in the right position.

**Picture 2: the flexible IUD eases through the cervix.**

**Picture 3: the IUD is inserted in the womb.**

**Picture 4: the insertion tube is removed.**

As soon as the Gynetics® IUD Cu250/ Cu375 is inserted, you are protected against pregnancy.

**Are there any side -effects?**

Shortly after insertion of the Gynetics® IUD Cu250/ Cu375, you can experience pain or cramps, which usually only lasts a short while. The first three menstruation periods might differ from your normal ones. They might be more painful and longer and sometimes there may be some bleeding in the periods between menstruation. After two or three months, it should be normal again. In case of abnormal discharge, pain during intercourse, abdominal cramps or fever, you should see your physician to exclude the possibility of pelvic inflammatory disease (PID).

Despite the use of an IUD, there is a slight chance of pregnancy. If you are over two weeks late, contact your physician. There might also be a slightly greater risk of an extra-uterine pregnancy.

Two months after insertion, your physician checks if the IUD is still properly in place. It is important to go for regular check-ups.

The IUD should be replaced every 3 years (Cu250) or 5 years (Cu375).

Your physician has explained to you how to check the presence of the nylon threads in the vagina. You should check the threads regularly with clean hands.

**What is the effect of the Gynetics® IUD on fertility?**

If you intend to get pregnant, ask your physician to remove the Gynetics® IUD Cu250/ Cu375. After removal, pregnancy is possible.

**The sterile IUD is for single use only. Do not re-use the IUD.**

**Gynetics® IUDs are available on prescription.**


**Insertion should be performed by a physician.**




If you are wearing a copper IUD, you must inform the attending physician of this before undergoing an MRI scan.

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