GYNÉTICS®

SOF-T® Cu375

0120

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SOF-T® Cu375

Information on the intrauterine contraceptive device, a method to prevent pregnancy:

World-wide, millions of women use an intrauterine contraceptive device as a method to prevent pregnancy. There are different kinds of intrauterine contraceptive devices- or IUDs. IUD is the abbreviation for intra-uterine device. Studies show that the IUD is a reliable and safe method of anti-conception. The Gynétics® SOF-T® IUD is made of plastic and is partially wrapped with copper. The IUD is shaped according to the anatomy of the uterus.

Picture C: position of the SOF-T® IUD in the uterus

- (1) womb
- (2) neck of the womb
- (3) cervix
- (4) vagina
- (a) horizontal arms with occlusion pearls
- (b) stem with copper
- (c) nylon threads

How does the Gynétics® IUD SOF-T® Cu375 work?

Compared to other temporary and reversible methods, the IUD is almost as safe as the birth control pill. Although no method to prevent pregnancy is a 100% safe, a SOF-T® Cu375 IUD is very reliable when it comes to protecting you against pregnancy. An IUD is does not protect you from sexually transmitted diseases.

The mechanism of the IUD is still not completely clear. It is likely that the presence of the IUD disturbs the normal function of the male sperm cells, which prevents fertilisation of the egg. It is also likely that the copper effects the development of the egg, which prevents fertilisation of the egg. For this reason IUDs are not thought to induce abortion.

How is the Gynétics® IUD SOF-T® Cu375 inserted?

The best time for insertion of the Gynétics® SOF-T® IUD is during menstruation. Most women tolerate insertion well. Before your physician inserts the IUD, the womb is examined and the length ascertained. The physician has also examined you carefully to see whether the IUD is a suitable method to prevent pregnancy for you. The IUD is carefully inserted in the right position:

- Picture A: the IUD eases through the cervix.
- Picture B: the plunger pushes the SOF-T® into its proper place and the IUD is inserted in the womb. The insertion tube and plunger are removed and the nylon thread is cut at a length of 1-2 cm.

As soon as the Gynétics® SOF-T® IUD is inserted, you are protected against pregnancy.

Are there any side -effects?

Shortly after insertion of the SOF-T® Cu 375, you can experience pain or cramps, which usually only lasts a short while. The first three menstruation periods could differ from your normal ones. They might be more painful and longer and sometimes there may be some bleeding in the periods between menstruation. After two or three months it should be normal again. In case of abnormal discharge, pain during intercourse, abdominal cramps or fever, you should see your physician to exclude the possibility of pelvic inflammatory diseases (PID).

Despite the use of an IUD, there is a slight chance of pregnancy. If you are over two weeks late, contact your physician. There might also be a slightly greater risk of an extra-uterine pregnancy. Women with a SOF-T® Cu375 IUD are less likely to have an ectopic pregnancy; properly inserted, the occlusion pearls block the fallopian tubes at least partially. Two months after insertion, your physician checks if the IUD is still properly in place. It is important to go for regular check-ups. The IUD should be replaced every 5 years.

Your physician has explained to you how to check the presence of the nylon threads in the vagina. You should check the threads regularly with clean hands.

What is the effect of the Gynétics® SOF-T® IUD on fertility?

If you intend to get pregnant, ask your physician to remove the SOF-T® Cu375 IUD. After removal, pregnancy is possible.

The sterile IUD is for single use only. Do not re-use the IUD.

The Gynétics® SOF-T® IUD is available on prescription. Insertion should be performed by a physician


If you are wearing a copper IUD, you must inform the attending physician of this before undergoing an MRI scan.


Manufacturer:
Gynétics Medical Products n.v., Rembert Dodoenstraat 51, 3920 Lommel, Belgium

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