
#5164 ECHOSAMPLER

Saline Infusion Sonography Canula and Endometrial Curette

MEDICAL
DEVICE

1639

Important product information, read carefully before use

Disposable – discard after single use**Sterile: unless pouch is damaged or opened****Re-use may cause sexually transmitted infections****Caution: federal (USA) law restricts this device to sale by or on the order of a physician****INDICATIONS FOR USE**

The EchoSampler™ is a single-use, sterile, disposable curette for

- Detection of abnormalities of the uterus and endometrium using real-time sonography during injection of sterile saline into the uterine cavity
- obtaining histological biopsy of the glandular epithelium and superficial chorionic layers of the uterine endometrial wall or sample extraction of uterine menstrual content for any of the following:
 - Routine screening for early detection of endometrial carcinoma or other precancerous conditions which could make estrogen therapy inadvisable.
 - Evaluation of endometrial tissue response in patients receiving estrogen replacement therapy for menopausal symptoms.
 - Endometrial dating and evaluation of uterine pathology associated with infertility, luteal insufficiency, or functional metrorrhagia.
 - Identification of specific uterine pathogens by bacterial culturing of uterine samples.
 - Inadequate imaging of the endometrium by endovaginal sonography
 - Further evaluation of suspected abnormalities as seen on endovaginal sonography, including focal or diffuse endometrial thickening or debris
 - Preoperative and postoperative evaluation of the uterine cavity, especially with regard to uterine polyps, myomas and cysts
 - Congenital abnormalities and/or anatomic variants of the uterine cavity
 - Infertility and habitual abortion

CONTRAINDICATIONS

The EchoSampler™ should not be used when one or more of the following conditions exist:

- Pregnancy or suspicion of pregnancy.
- Acute pelvic inflammatory disease (PID) or recent treatment for PID.
- Untreated acute cervicitis, chronic cervicitis, or vaginitis, including bacterial vaginosis, until infection is controlled.

Pelvic organ tenderness should be assessed during the preliminary endovaginal sonogram. Active vaginal bleeding is not a contraindication to the procedure but may take the procedure more challenging.

**WARNINGS**

The EchoSampler™ should not be used to obtain an endometrial biopsy in patients with amenorrhoea unless a laboratory test has confirmed the absence of detectable circulating HCG levels.

The EchoSampler should be used by healthcare professionals only.

PRECAUTIONS

Prior to insertion of the EchoSampler™, the uterus should be carefully sounded to determine the degree of patency of the endocervical canal and the internal os, and the direction and depth of the uterine cavity. If resistance is felt at any time in the procedure, the EchoSampler™ should never be forced.

Slight lubrication of the sheath with a water soluble gel may aid the insertion.

The physician may elect to prescribe prophylactic antibiotics if the

patient routinely take these for other invasive procedures. If painful, dilated and/or obstructed fallopian tubes are found prior to saline infusion, and the patient is not taking prophylactic antibiotics, the examination should be delayed until treatment can be administered. In the presence of nontender hydrosalpinges, consideration may be given to administer antibiotics at the time of the examination.

ADVERSE REACTIONS

The adverse reactions cited as having been reported to occur are not listed in any order of frequency or severity. Reported adverse reactions from endometrial sampling procedures include perforation of the uterus, pain and cramping, uterine spasm, vasovagal syncope, and vaginal bleeding.

**Manufactured by:**Gynetics Medical Products N.V.
Rembert Dodoensstraat 51
3920 Lommel, Belgium**HOW SUPPLIED**

#5164 Echosampler

PACKAGED

5 boxes of 25 individually sterile packed, disposable devices.

DIRECTIONS FOR USE

- Perform a preliminary unenhanced endovaginal sonography with measurements of endometrium and evaluation of the uterus and ovaries
- Prepare the vagina and cervix using currently accepted antiseptic techniques for intrauterine procedures, use of aseptic technique during the entire procedure is essential.
- With an open-sided speculum in place, gently insert a sterile uterine sound to determine the depth and direction of the uterine canal. If the uterus is anteverted or retroverted, it may be advisable to use very fine forceps or a tenaculum to correct the angulation and stabilize the cervix.
- After sounding the uterus, the depth and direction of the uterine canal should be noted. With the saline-filled syringe attached to the Luer Lock connection the EchoSampler™ is gently inserted through the cervical canal into the uterine cavity. If resistance is encountered, no attempt should be made to force the insertion. In patients with an extremely dry or narrow cervical canal, slight lubrication of the sheath with a water soluble gel may aid the insertion.
- Remove the open-sided speculum and tenaculum while keeping the EchoSampler™ in position.
- Insert a high-frequency vaginal transducer covered by a protective sheath with sterile coupling gel caudal to the EchoSampler™ and visualise the uterus
- With the sheath correctly positioned in the uterine cavity, the piston of the syringe should be pressed administering the sterile saline into the uterine cavity under real-time sonographic imaging. The echogenic sheath of the EchoSampler™ will indicate the position of the tip in relation to the uterine cavity. Representative images with a complete survey of the uterine cavity are obtained as necessary for diagnostic evaluation
- Pulling back the piston of the syringe for about 10 cc will create the maximum negative pressure within the sheath of the EchoSampler™ and result in an optimal tissue sample, the sheath should be rolled between the fingers while simultaneously moving it laterally as well as back and forth inside the uterus three (3) to four (4) times for comprehensive sampling.
- The EchoSampler™ should be removed gently from the patient.
- The distal tip of the sheath should be examined for the presence of a uterine mucosa sample.
- To expel the sample into the transport container, the distal tip of the EchoSampler™ is cut just below the side hole and the piston is pushed into the syringe to expel the sample from the sheath.

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