#2000 SEMTRAC SET

Soft Embryo Transfer Catheter

HOW SUPPLIED:

REF #2000 Semtrac Set 5	REF #2000 Set 5 TIS
REF #2000 Semtrac Set 5 SP	REF #2000 Set 5 TIS SP
REF #2000 Semtrac Set 4.2	REF #2000 Set 4.2 TIS
REF #2000 Semtrac Set 4.2 SP	REF #2000 Set 4.2 TIS SP
REF #2000 Semtrac Set 5 A	REF #2000 Semtrac Set 5 AB
REF #2000 Semtrac Set 5 Trial	REF #2000 Semtrac Set C5
REF #2000 B Set 5	REF #2000 Semtrac Set Short
REF #2000 B Set 4.2	

Important product information, read carefully before use

Disposable – discard after single use

Re-use may cause sexually transmitted infections

Caution: federal (USA) law restricts this device to sale by or on the order of a physician.

The trial devices are intended to support the use of an approved embryo transfer catheter. These tools are not intended to directly transfer gametes/embryos into the uterine cavity; rather, they are to aid in traversing a difficult or stenotic cervix or to perform a trial/mock transfer to better identify the uterine environment.

Sterile: unless pouch is damaged or opened

DESCRIPTION:

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers. It is a three piece catheter with markings to indicate depth and has a handle with luer lock. Available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance.

INDICATIONS FOR USE:

Introduction of embryos into the uterine cavity.

CONTRAINDICATIONS:

Do not use in the presence of cervical infection.

Do not use in the presence or after recent Pelvic Inflammatory disease.

Do not use for intrafallopian procedures.

**WARNINGS:**

Not sterile if pouch is damaged or opened, discard immediately.

The #2000 SEMTRAC SET should be used by healthcare professionals only.

Forcing the catheter may result in damage to the endometrial tissue and bleeding.

PRECAUTIONS:

By virtue of its flexibility and small diameter, the #2000 Semtrac Set can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroflexion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

In no case should the device be forced against digitally felt resistance. If passage of the #2000 Semtrac Set through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered.

ADVERSE REACTIONS:

No adverse reactions associated with the use of the #2000 Semtrac Set have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping and superficial haemorrhage may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

1639

**Manufactured by:**

Gynetics Medical Products N.V.
Rembert Dodoensstraat 51
3920 Lommel, Belgium

PACKAGED:

#2000 Semtrac Set 5/ #2000 Semtrac Set 4.2/#2000 Set 5 TIS/#2000 Set 4.2 TIS:

3 boxes of 25 individually sterile packed, disposable devices and 5 mandrels.

#2000 Semtrac Set 5 Trial:

3 boxes of 25 individually sterile packed, disposable devices

#2000 Semtrac Set 5 SP/ #2000 Semtrac Set 4.2 SP/#2000 Set 5 TIS SP/ #2000 Set 4.2 TIS SP:

3 boxes of 25 individually sterile packed, disposable sets of devices and 5 mandrels. Packaged in a double pouch.

#2000 Semtrac Set 5 A:

3 boxes of 25 individually sterile packed guiding catheters.

#2000 Semtrac Set Short:

3 boxes of 25 sterile packed sets of guiding and loading catheter

#2000 B Set 4.2/ #2000 B Set 5:

5 mandrels per 25 #2000 Semtrac Set 5/ #2000 Semtrac Set 4.2

#2000 Semtrac Set 5 AB:

5 individually sterile packed sets per 25 #2000 Semtrac Set C5

DIRECTIONS FOR USE FOR THE #2000 SEMTRAC SET 5 (SP) /#2000 SET 5 TIS (SP)/ #2000 SEMTRAC SET 4.2 (SP) / #2000 SET 4.2 TIS (SP)/ #2000 SEMTRAC SET 5 A/#2000 SEMTRAC SET C5/ #2000 SEMTRAC SET SHORT/ #2000 SEMTRAC SET 5 AB:

With the patient in the dorsal lithotomy position expose the uterine cervix to view with a vaginal speculum.

1. Having previously determined the depth and direction of the patient's uterus, gently advance the guiding catheter (white handle), into the external cervical os and advance it near the cervical canal into the uterine cavity until resistance is felt. Now retract the catheter approximately 1.0 cm by following the guides. Gently rotate the catheter to remove any excess mucus from the tip.

Note: If necessary, the use of the mandrel (blue handle) can facilitate insertion of the catheter. This by preforming and reinforcing the guiding catheter.

2. Remove the Mandrel from the guiding catheter (if applicable)
3. Load the embryo(s) into the loading catheter (pink handle), following your standard protocol.
4. Insert loading catheter (pink handle) into the guiding catheter (white handle).
5. The markings on the proximal end of the loading catheter (pink handle) indicate the distance the loading cannula has advanced beyond the tip of guiding catheter (white handle).
6. The distal tip of the catheter should then be at the uterine fundus. Note: The tip of the TIS loading catheters can be viewed under ultrasound.
7. Note: If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. In no case should force be used against resistance to achieve insertion.
8. Slowly release the embryos into the uterus.
9. Following the complete expulsion of the embryos, slowly withdraw the catheter from the uterus and discard it. Remove the speculum and allow the patient to rest.
10. Once the set is removed from the patient, inspect the loading catheter (pink handle) to verify embryo release.
11. Discard the catheter in accordance with local regulation.

DIRECTIONS FOR USE FOR THE #2000 SEMTRAC SET 5**TRIAL:**

1. Place patient in lithotomy or knee/chest position.
2. Insert vaginal speculum to expose cervix.
3. Use malleable/flexible obturator with guiding catheter to traverse cervix if necessary. Remove malleable/flexible obturator prior to placement of trial/mock catheter. (May be performed under ultrasound guidance if desired.)
4. Introduce trial/mock catheter through guiding catheter.
5. Remove trial/mock catheter from guiding catheter.
6. If so desired, move on to embryo transfer using approved catheter and refer to instructions for use for that approved catheter.
7. Upon completion of trial/mock procedure, remove the catheter/catheters, and discard.


Catalogue number

Batch code

Manufacturer

Temperature limit

Consult instructions for use

CE marking of conformity

Use-by date

Sterilization using ethylene oxide

Do not use if package is damaged

Do not reuse

Medical device

Do not resterilize

#2000 SEMTRAC SET

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REF #2000 Semtrac Set 5	REF #2000 Set 5 TIS
REF #2000 Semtrac Set 5 SP	REF #2000 Set 5 TIS SP
REF #2000 Semtrac Set 4.2	REF #2000 Set 4.2 TIS
REF #2000 Semtrac Set 4.2 SP	REF #2000 Set 4.2 TIS SP
REF #2000 Semtrac Set 5 A	REF #2000 Semtrac Set 5 AB
REF #2000 Semtrac Set 5 Trial	REF #2000 Semtrac Set C5
REF #2000 B Set 5	REF #2000 Semtrac Set Short
REF #2000 B Set 4.2	

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POKYNY PRO POUŽITÍ SAD #2000 SEMTRAC SET 5 (SP) /#2000 SET 5 TIS (SP) / #2000 SEMTRAC SET 4.2 (SP) / #2000 SET 4.2 TIS (SP) / #2000 SEMTRAC SET 5 A/#2000 SEMTRAC SET C5/ #2000 SEMTRAC SET SHORT/ #2000 SEMTRAC SET 5 AB:

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Blødt embryo-transferkateter

LEVERES SÅDAN:

REF #2000 Semtrac Set 5	REF #2000 Set 5 TIS
REF #2000 Semtrac Set 5 SP	REF #2000 Set 5 TIS SP
REF #2000 Semtrac Set 4.2	REF #2000 Set 4.2 TIS
REF #2000 Semtrac Set 4.2 SP	REF #2000 Set 4.2 TIS SP
REF #2000 Semtrac Set 5 A	REF #2000 Set 4.2 TIS SP
REF #2000 Semtrac Set 5 Trial	REF #2000 Semtrac Set 5 AB
REF #2000 B Set 5	REF #2000 Semtrac Set C5
REF #2000 B Set 4.2	REF #2000 Semtrac Set Short

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