#4000 THE TULIP

Embryo Transfer catheter

HOW SUPPLIED:

REF #4000 The Tulip	REF #4000 S Memo
REF #4000 The Tulip DP	REF #4000 L Memo
REF #4000 The Tulip Memo	REF #4000 LS Memo
REF #4000 The Tulip Memo A	REF #4000 S Memo TIS
REF #4000 The Tulip A	REF #4000 L Memo TIS
REF #4000 Tulip S	REF #4000 LS Memo TIS
REF #4000 Tulip L	REF #4000 Memo B L-LS
REF #4000 Tulip LS	REF #4000 The Tulip Memo B
REF #4000 Tulip S TIS	REF #4000 Tulip LS TIS
REF #4000 Tulip L TIS	

Important information, read carefully before use

Disposable – discard after single use

Re-use may cause sexually transmitted infections

Caution: federal (USA) law restricts this device to sale by or on the order of a physician

Sterile: unless pouch is damaged or opened

DESCRIPTION:

The Tulip is a reshaped embryo transfer catheter with a tulip shaped bulb for difficult cervical passage, consisting of two parts: a guiding catheter (A) and a loading catheter (C).

INDICATIONS FOR USE:

Introduction of embryos into the uterine cavity.

CONTRAINDICATIONS:

Do not use in the presence of cervical infection.

Do not use in the presence or after recent Pelvic Inflammatory disease.

Do not use for intrafallopian procedures.

WARNINGS:

Not sterile if pouch is damaged or opened, discard immediately.

The Tulip should be used by healthcare professionals only

In no case should the device be forced against digitally felt resistance. If passage of #4000 The Tulip / #4000 The Tulip Memo through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered.

PRECAUTIONS:

By virtue of its flexibility and small diameter, #4000 The Tulip / #4000 The Tulip Memo can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine ante flexion or retro flexion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

ADVERSE REACTIONS:

No adverse reactions associated with the use of #4000 The Tulip / #4000 The Tulip Memo have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping and superficial haemorrhage may be expected. In every case the direction for use should be followed; taking note of contraindications, warnings and precautions.

1639

Manufactured by:

Gynétics Medical Products N.V.
Rembert Dodoensstraat 51
3920 Lommel, Belgium

PACKAGED:

- #4000 The Tulip / #4000 The Tulip Memo / #4000 Tulip L (TIS) / #4000 Tulip S (TIS) / #4000 Tulip LS (TIS) / #4000 S Memo (TIS) / #4000 L Memo (TIS) / #4000 LS Memo (TIS) :
3 boxes of 25 individually sterile packed, disposable sets of devices.

- #4000 The Tulip Memo B:
5 mandrels per 25 #4000 The Tulip Memo
- #4000 Memo B L-LS:
5 mandrels per 25 #4000 S Memo / #4000 L Memo / #4000 LS Memo / #4000 S Memo TIS / #4000 L Memo TIS / #4000 LS memo TIS.
- #4000 The Tulip A:
3 boxes of 25 individually sterile packed guiding catheters
- #4000 The Tulip Memo A:
3 x 25 individually sterile packed guiding catheters
- #4000 Tulip DP:
3 x 25 Sterile packed sets of a individually packed guiding catheter and loading catheter

DIRECTIONS FOR USE:

1. With the patient in the dorsal lithotomy position expose the uterine cervix to view with a vaginal speculum.
 2. Having previously determined the depth, set positioner on the guiding catheter (white handle).
 3. Having previously determined direction of the patient's uterus, put the mandrel, (blue handle), in guiding catheter (white handle) to pre-shape The Tulip A / The Tulip Memo A.
 4. Gently advance the catheter into the external cervical os and advance it through the cervical canal near the uterine cavity until resistance is felt.
 5. Remove the mandrel (blue handle) when position is reached
 6. Load the embryo(s) into the loading catheter (pink handle) following your standard protocol. Insert the loading catheter. The marking on the proximal end of the tube indicates when the distal tip of the guiding catheter is reached.
 7. The distal tip of the loading catheter should then be at the uterine fundus. **Note:** The tip of the TIS loading catheters can be viewed under ultrasound.
- Note:** If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. In no case should force be used against resistance to achieve insertion.
8. Slowly release the embryos into the uterus.
 9. Following the complete expulsion of the embryos, slowly withdraw the catheter from the uterus and discard it. Remove the speculum and allow the patient to rest.
 10. Once the set is removed from the patient, inspect The Tulip C / #4000 The Tulip Memo C to verify embryo release.
 11. Safely discard after use in accordance with local regulation.

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