#4551 FAS

Follicle aspiration needle

HOW SUPPLIED:**REF** #4551 FAS

Important product information, read carefully before use.

Disposable – discard after single use.**Re-use may cause sexually transmitted infections.****Caution: federal (USA) law restricts this device to be used by a physician.****Sterile: unless pouch is damaged or opened.****DESCRIPTION:**

The device is indicated for use in ultrasound guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.

The device is available with:

- a single lumen (without extra tap for flushing)
- a single lumen flush (with extra tap for flushing)
- a double lumen (needle within needle, each with a separate tap, for flushing and suction)

INDICATIONS FOR USE:

Transvaginal harvesting of oocytes from follicles for IVF and ICSI.

CONTRAINDICATIONS:

Do not use in the presence of cervical infection, after recent pelvic inflammatory disease, sexually transmitted disease or who is currently pregnant. Do not use if the package is opened/damaged; discard immediately.

**WARNINGS:**

Please refer to published medical literature before performing this procedure. The device should be used by healthcare professionals only. The device should not be used in case of resistance to achieve insertion.

ADVERSE REACTIONS:

No adverse reactions associated with the use of the device have been reported. However, as with any invasive instrument cramping and mild bleeding may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

**Manufactured by:**

Gynétics Medical Products N.V.
Rembert Doddeensstraat 51
3920 Lommel, Belgium

PACKAGED:

2 boxes of 15/20/25 individually sterile packed, disposable devices.

DIRECTIONS FOR USE:

1. Prepare the patient for oocyte aspiration according to standard clinical practice
2. Check the expiry date and connectivity of all parts.
3. Connect the needle to the pump or syringe – depending on the type of needle and insert the (silicone) cork into a (sterile) vial if applicable
4. Insert the needle into the ovarian follicle following your standard protocol by using ultrasound guidance to locate the follicles.
5. Perform follicle aspiration with vacuum pump or syringe. The pump pressure used with a specific needle type and size is at the discretion of the clinician performing the procedure.
If flushing is necessary: proceed with one of the following:
 - a) Single Lumen Flush:
Attach a syringe with sterile media to the luer compatible handle and flush.
 - b) Double Lumen:
Attach a syringe with sterile media to the luer compatible handle and flush. Continue to aspirate the follicle until it is empty
6. Remove the needle, repeat steps 4 and 5 on the remaining follicles
7. After aspirating the desired number of oocytes, withdraw the needle from the uterus.
8. Allow the patient to rest.
9. Discard the needle in accordance with local regulation.

Catalogue number

Batch code

Manufacturer

Temperature limit

Consult instructions for use

CE marking of conformity

Use-by date

Sterilized using irradiation

Do not use if package is damaged

Do not reuse

Medical device

Do not resterilize

Single sterile barrier system

Single sterile barrier system with protective packaging

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