
**#2000 SEMTRAC SET**

Soft Embryo Transfer Catheter

**HOW SUPPLIED:**

	#2000 Semtrac Set 5		#2000 Set 5 TIS
	#2000 Semtrac Set 5 SP		#2000 Set 5 TIS SP
	#2000 Semtrac Set 4.2		#2000 Set 4.2 TIS
	#2000 Semtrac Set 4.2 SP		#2000 Set 4.2 TIS SP
	#2000 B Set 5		
	#2000 B Set 4.2		

Important product information, read carefully before use

**Disposable – discard after single use**

**Re-use may cause sexually transmitted infections**

**Caution: federal (USA) law restricts this device to sale by or on the order of a physician.**

**Sterile: unless pouch is damaged or opened.**

**DESCRIPTION:**

The device is a soft embryo transfer catheter set for difficult embryo transfers. It is a two or three piece catheter with markings to indicate depth and has a handle with luer lock. The catheter is available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance.

**INDICATIONS FOR USE:**

Introduction of embryos into the uterine cavity.

**CONTRAINDICATIONS:**

Do not use in the presence of cervical infection, after recent pelvic inflammatory disease, sexually transmitted disease or who is currently pregnant. Do not use if the package is opened/damaged; discard immediately.

**WARNINGS:**

Please refer to published medical literature before performing this procedure. The device should be used by healthcare professionals only. Forcing the catheter may result in damage to the endometrial tissue and bleeding.

**PRECAUTIONS:**

By virtue of its flexibility and small diameter, the device can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroflexion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized). In no case should the device be forced against digitally felt resistance. If passage of the device through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered.

**ADVERSE REACTIONS:**

No adverse reactions associated with the use of the device have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping and superficial haemorrhage may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

1639

**Manufactured by:**

Gynetics Medical Products N.V.  
Rembert Dodoensstraat 51  
3920 Lommel, Belgium

**PACKAGED:**

**#2000 Semtrac Set 5/ #2000 Semtrac Set 4.2/ #2000 Set 5 TIS/ #2000 Set 4.2 TIS:**

3 boxes of 25 individually sterile packed, disposable devices and 5 mandrels.

**#2000 Semtrac Set 5 SP/ #2000 Semtrac Set 4.2 SP/ #2000 Set 5 TIS SP/ #2000 Set 4.2 TIS SP:**

3 boxes of 25 individually sterile packed, disposable sets of devices.

**#2000 B Set 4.2/ #2000 B Set 5:**

Mandrel individually sterile packed

**DIRECTIONS FOR USE:**

1. Check the 'expiry date' of the device
2. Position the patient in the lithotomy position. Insert vaginal speculum to expose cervix. Determine the depth and direction of the patient's uterus.
3. Gently advance the guiding catheter (white handle), into the external cervical os and advance it near the cervical canal into the uterine cavity until resistance is felt. Gently rotate the catheter to remove any excess mucus from the tip.

**Note:** If necessary, the use of the mandrel (blue handle) can facilitate insertion of the catheter. This by preforming and reinforcing the guiding catheter.

4. Remove the Mandrel from the guiding catheter (if applicable)
5. Attach the syringe to the luer lock of the loading catheter (pink handle). Load the embryo(s) into the loading catheter, following your standard protocol.
6. Insert loading catheter into the guiding catheter.
7. The markings on the proximal end of the loading catheter indicate the distance of the loading cannula that has advanced beyond the tip of the guiding catheter.
8. The distal tip of the catheter should then be at the uterine fundus.

**Note:** The tip of the TIS loading catheters can be viewed under ultrasound.

**Note:** If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. In no case should force be used against resistance to achieve insertion.

9. Slowly press the plunger of the syringe and expel the embryos into the uterine cavity.
10. Following the complete expulsion of the embryos, slowly withdraw the catheter(s) from the uterus. Inspect the loading catheter (pink handle) to verify embryo release.
11. Remove the speculum and allow the patient to rest.
12. Discard the catheter in accordance with local regulation.

Catalogue number

Batch code

Manufacturer

Temperature limit

Consult instructions for use

CE marking of conformity

Use-by date

Sterilization using ethylene oxide

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