#2000 SEMTRAC C

Soft Embryo Transfer Catheter

HOW SUPPLIED:

-  #2000 Semtrac C
-  #2000 C TIS
-  #2000 Semtrac Set C5
-  #2000 Semtrac Set 5 AB

Important product information, read carefully before use.

Disposable – discard after single use. Re-use may cause sexually transmitted infections. Caution: federal (USA) law restricts this device to be used by a physician. Sterile: unless pouch is damaged or opened

DESCRIPTION:

The device is a soft embryo transfer catheter. It is a one piece catheter with markings to indicate depth and has a handle with luer lock. #2000 Semtrac C is available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance. Intended patient population: mature women.

INDICATIONS FOR USE:

Introduction of embryos into the uterine cavity.

CONTRAINDICATIONS:

Do not use in the presence of cervical infection, after recent pelvic inflammatory disease, sexually transmitted disease or who is currently pregnant. Do not use if the package is opened/damaged; discard immediately.

WARNINGS:

Please refer to published medical literature before performing this procedure. The device should be used by healthcare professionals only. Forcing the catheter may result in damage to the endometrial tissue and bleeding. If passage of the device through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS:

By virtue of its flexibility and small diameter, the device can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroflexion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

ADVERSE REACTIONS:

No adverse reactions associated with the use of the device have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

1639

Manufactured by:
Gynetics Medical Products N.V.
Rembert Dodoensstraat 51
3920 Lommel, Belgium

PACKAGED:

#2000 Semtrac C/ #2000 C TIS

3 boxes of 25 individually sterile packed, disposable loading catheters

#2000 Semtrac Set C5 + 5 per 25 #2000 Semtrac Set 5 AB

3 boxes of 25 individually sterile packed, disposable sets of devices

DIRECTIONS FOR USE:

1. Check the expiry date of the device

2. Position the patient in the lithotomy position. Insert vaginal speculum to expose cervix. Determine the depth and direction of the patient's uterus,
3. Attach the syringe to the luer lock of the loading catheter (pink handle). Load the embryo(s) into the loading catheter, following your standard protocol.
4. Gently advance the loading catheter, into the external cervical os and advance it through the cervical canal into the uterine cavity until resistance is felt. The distal tip of the catheter should then be at the uterine fundus. Gently rotate the catheter to remove any excess mucus from the tip.

Note: The tip of the TIS loading catheters can be viewed under ultrasound.

Note: Use #2000 Semtrac 5AB as assistance for #2000 Semtrac set C5. If necessary, the use of the mandrel (blue handle) can facilitate insertion of the guiding catheter. This by preforming and reinforcing the guiding catheter.

Note: If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. The device should not be used in case of resistance to achieve insertion.

5. Slowly press the plunger of the syringe and expel the embryos into the uterine cavity.
6. Following the complete expulsion of the embryos, slowly withdraw the catheter(s) from the uterus. Inspect the loading catheter (pink handle) to verify embryo release.
7. Remove the speculum and allow the patient to rest.
8. Discard the catheter in accordance with local regulation.

 Catalogue number

 Batch code

 Manufacturer

 Temperature limit
0°C / 32°F to 30°C / 86°F

 Consult instructions for use

 CE marking of conformity

 Use-by date

 Sterilization using ethylene oxide

 Do not use if package is damaged

 Do not reuse

 Medical device

 Do not resterilize

 Single sterile barrier system

 Single sterile barrier system with protective packaging

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