**#5165-LN970**

Gel Instillation Sonohysterography Catheter

**HOW SUPPLIED:**

**REF** #5165-LN970

Important product information, read carefully before use

**Disposable – discard after single use. Re-use may cause sexually transmitted infections. Caution: federal (USA) law restricts this device to sale by or on the order of a physician. Sterile: unless pouch is damaged or opened**

**DESCRIPTION:**

The Sterile GIS catheter: gel instillation sonohysterography catheter is assembled of tubing, an inject moulded handle and a inject moulded accessory. Intended patient population: women.

**INDICATIONS FOR USE:**

Introduction of a gel or foam into the uterine cavity and/or fallopian tubes for all indications for sonohysterography.

The GIS Catheter is intended to be used by Health care professionals in the field of gynecology and fertility

**CONTRAINDICATIONS:**

Do not use in the presence of active pelvic infection, sexually transmitted disease, profuse bleeding or pregnancy. Nor in case of allergy to any of the constituents of ExEm® Gel (hydroxyethylcellulose, glycerol and purified water).

**PRECAUTIONS AND WARNINGS:**

The GIS Catheter should be considered as contaminated medical waste after use. Do not use if package of GIS-catheter has been previously opened or damaged. The GIS-Catheter is intended for single use only. Do not reuse, re-sterilize or autoclave. Reuse of single use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. The manufacturer will not be responsible for any direct, incidental or consequential damages resulting from re-sterilization or reuse. Do not use the GIS Catheter after the expiry date. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**ADVERSE REACTIONS:**

No adverse reactions associated with the use of the GIS catheter have been reported. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

1639

**Manufactured by:**  
Gynetics Medcal Products N.V.,  
Rembert Dodoensstraat 51,  
3920 Lommel, Belgium

**PACKAGED:**

The individually sterile packed, disposable #5165-LN970 catheters are packed in a box.

**DIRECTIONS FOR USE FOR EXEM® FOAM:**

1. Connect a syringe, filled with a foam suitable for sonohysterosalpingography, to the purple handle of the inner tubing/catheter.

2. Fill the cervical applicator with foam prior to introduction
3. Insert the cervical applicator of the catheter into the cervical canal of the uterus until a certain resistance is felt.
4. Slowly infuse small amounts of foam and complete the ultrasound procedure.
5. Remove the catheter and cervical applicator.
6. Discard as medical device in accordance with local regulation


**DIRECTIONS FOR USE FOR GIS-KIT**

1. Connect a syringe filled with ExEm® Gel suitable for sonohysterography, to the purple handle of the inner tubing/catheter (Fig. 1).
2. Fill the cervical applicator with ExEm® Gel until the tubing and applicator are free of air bubbles (Fig. 2).
3. Insert the cervical applicator into the cervical canal of the uterus until a certain resistance is felt and instill the ExEm® Gel into the uterus until the patient experiences a certain feeling of pain (Fig. 3).
4. Disconnect the syringe and leave the cervical applicator in place by pulling the syringe (and connected inner tubing/catheter) and holding the outer tubing of the catheter; some drops of ExEm® Gel may leak (Fig. 4).
5. Complete the ultrasound procedure.
6. Remove the cervical applicator by pulling the cord, this can be done by the ultrasonographer or the patient (Fig. 5).
7. Discard as medical device in accordance with local regulation

**REF** Catalogue number

**LOT** Batch code

 Manufacturer

 Temperature limit  
30°C / 86°F  
0°C / 32°F

 Consult instructions for use

 CE marking of conformity

 Use-by date

 Sterilization using ethylene oxide

 Do not use if package is damaged

 Do not reuse

 Medical device

 Do not resterilize


 Single sterile barrier system

 Single sterile barrier system with protective packaging

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1639

**Manufactured by:**  
Gynetics Medical Products N.V.,  
Rembert Dodoensstraat 51,  
3920 Lommel, Belgium

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