**#4000 THE TULIP**

Embryo Transfer catheter

**HOW SUPPLIED:**

 #4000 The Tulip	 #4000 S Memo
 #4000 The Tulip Memo	 #4000 L Memo
 #4000 Tulip S	 #4000 LS Memo
 #4000 Tulip L	 #4000 S Memo TIS
 #4000 Tulip LS	 #4000 L Memo TIS
 #4000 Tulip S TIS	 #4000 LS Memo TIS
 #4000 Tulip L TIS	 #4000 Memo B L-LS
 #4000 Tulip LS TIS	 #4000 The Tulip Memo B

Important information, read carefully before use.

**Disposable – discard after single use. Re-use may cause sexually transmitted infections. Caution: federal (USA) law restricts this device to sale by or on the order of a physician. Sterile: unless pouch is damaged or opened.**

**DESCRIPTION:**

The device is a (preshaped) embryo transfer catheter with a tulip shaped bulb for difficult cervical passage. It is a two or three piece catheter with markings to indicate depth, has a handle with luer lock and a positioner for depth adjustment. The catheter is available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance. Intended patient population: mature women.

**INDICATIONS FOR USE:**

Introduction of embryos into the uterine cavity.

**CONTRAINDICATIONS:**

Do not use in the presence of cervical infection, after recent pelvic inflammatory disease, sexually transmitted disease or who is currently pregnant. Do not use if the package is opened/damaged; discard immediately.

**WARNINGS:**

Please refer to published medical literature before performing this procedure. The device should be used by healthcare professionals only. Forcing the catheter may result in damage to the endometrial tissue and bleeding. If passage of the device through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**PRECAUTIONS:**

By virtue of its flexibility and small diameter, the device can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroversion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

**ADVERSE REACTIONS:**

No adverse reactions associated with the use of the device have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping and superficial haemorrhage may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

1639

**Manufactured by:**

Gynetics Medical Products N.V.  
Rembert Dodoensstraat 51  
3920 Lommel, Belgium

**PACKAGED:**

**#4000 The Tulip / #4000 Tulip L (TIS) / #4000 Tulip S (TIS) / #4000 Tulip LS (TIS):**

3 boxes of 25 individually sterile packed, disposable sets of devices.

**#4000 The Tulip Memo / #4000 S Memo (TIS) / #4000 L Memo (TIS) / #4000 LS Memo (TIS):**

3 boxes of 25 individually sterile packed, disposable devices and 5 mandrels.

**#4000 The Tulip Memo B/ #4000 Memo B L-LS:**

Mandrel individually sterile packed

**DIRECTIONS FOR USE:**

1. Check the expiry date of the device.
2. Position the patient in the lithotomy position. Insert vaginal speculum to expose the cervix. Determine the depth and direction of the patient's uterus.
3. Having previously determined the depth, set positioner on the guiding catheter (white handle).
4. Gently advance the guiding catheter into the external cervical os and advance it through the cervical canal near the uterine cavity until resistance is felt. Gently rotate the catheter to remove any excess mucus from the tip.

**Note:** If necessary, the use of the mandrel (blue handle) can facilitate insertion of the guiding catheter. This by preforming and reinforcing the guiding catheter.

5. Remove the mandrel (blue handle) when position is reached (if applicable)
6. Attach the syringe to the luer lock of the loading catheter (pink handle). Load the embryo(s) into the loading catheter following your standard protocol.
7. Insert the loading catheter in the guiding catheter.
8. The markings on the proximal end of the loading catheter indicate the distance of the loading cannula that has advanced beyond the tip of the guiding catheter.
9. The distal tip of the loading catheter should then be at the uterine fundus. **Note:** The tip of the TIS loading catheters can be viewed under ultrasound.

**Note:** If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. The device should not be used in case of resistance to achieve insertion.

10. Slowly press the plunger of the syringe and expel the embryos into the uterine cavity.
11. Following the complete expulsion of the embryos, slowly withdraw the catheter(s) from the uterus. Inspect the loading catheter (pink handle) to verify embryo release.
12. Remove the speculum and allow the patient to rest.
13. Discard the catheters in accordance with local regulation.

Catalogue number

Batch code

Manufacturer

Temperature limit

Consult instructions for use

CE marking of conformity

Use-by date

Sterilized using irradiation

Do not use if package is damaged

Do not reuse

Medical device

Do not resterilize

Single sterile barrier system

Single sterile barrier system with protective packaging

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