#4220 IUI/#4220 IUI MEMO/#4220 PLUS

Intra uterine Insemination Cannula

HOW SUPPLIED: #4220 IUI #4220 IUI Memo #4220 Plus

Important product information-Read carefully before use.

Disposable – discard after single use.**Re-use may cause sexually transmitted infections.****Caution: federal (USA) law restricts this device to be used by a physician.****Sterile: unless pouch is damaged or opened.****DESCRIPTION:**

The device is an intrauterine insemination cannula with two lateral distal ports, for #4220 IUI and #4220 IUI Memo for optimal sperm distribution. The #4220 IUI Memo has a memo function which proactively enables the distal tip to bend according to cervical flexion. The #4220 plus has an open tip. The #4220 IUI (Memo) is a one piece catheter. The #4220 IUI Plus is a two piece catheter which includes additionally a mandrel. Intended patient population: individuals, male or female with infertility disorder undergoing medical assisted reproduction procedures.

INDICATIONS FOR USE:

Intra uterine insemination

CONTRAINDICATIONS:

Do not use in the presence of cervical infection, after recent pelvic inflammatory disease, sexually transmitted disease or who is currently pregnant. Do not use if the package is opened/damaged; discard immediately.

**WARNINGS:**

Always use washed spermatozoa when performing intrauterine (including intracervical) insemination. Please refer to published medical literature for methods of preparing spermatozoa for intrauterine insemination before performing this procedure. If passage of the device through the cervical canal is not found to be easily achievable, the device should not be forced. The device should be used by healthcare professionals only trained in medical assisted reproduction procedures. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS:

By virtue of its flexibility and small diameter, the device can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroflexion present in individual patients. (This is preferably done during a previous examination and not on the day of the insemination when disturbance of the uterus should be minimized).

ADVERSE REACTIONS:

No adverse reactions associated with the use of the device have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

1639

**Manufactured by:**

Gynetics Medical Products N.V.
Rembert Dodoensstraat 51
3920 Lommel, Belgium

DIRECTIONS FOR USE:

1. Check the expiry date of the device.
2. Position the patient in the lithotomy position. Insert vaginal speculum to expose the cervix. Determine the depth and direction of the patient's uterus.
3. A) Firmly attach a syringe (preferably 5ml or 3ml) with the plunger fully depressed, to the handle at the proximal end of the device. Draw 1.0 ml of air into the syringe. Then place the distal tip of the device into the specimen of washed spermatozoa so that the side port is immersed in the specimen. Now aspirate the desired quantity of the specimen into the lumen of the catheter by retraction of the syringe plunger.
B) Load the syringe according to standard clinical practice. To inseminate attach the device at the tip of the syringe.
4. Adjust the positioner corresponding to the anatomical proportions determined before.

Note: If necessary, the use of the mandrel (blue handle) can facilitate insertion of the catheter. This by preforming and reinforcing the catheter.

5. Remove the mandrel (blue handle) when position is reached (if applicable)
6. Insert the catheter until the positioner is being on the outer uterine orifice. Rotate the catheter to bring the arrows on the handle to the side. This orients the two lateral ports at the distal end of the catheter within the uterine cavity toward the uterine ostia.

Note: The small OD of the distal segment of the device and the added control provided by its more rigid proximal segment usually permits its easy insertion. If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. The device should not be used in case of resistance to achieve insertion.

7. Inject the suspension with the spermatozoa slowly into the cave uteri.
8. Following the complete expulsion of the specimen, slowly withdraw the device from the uterus and discard it in accordance with local regulation
9. Remove the speculum and allow the patient to rest.

Catalogue number

Use-by date

Batch code


Sterilized using irradiation

Manufacturer

Do not reuse

Temperature limit

Do not use if package is damaged

Consult instructions for use

Medical device

CE marking of conformity

Do not resterilize

Single sterile barrier system

Single sterile barrier system with protective packaging

Unique device identifier

Keep dry

Date of manufacture

PACKAGED:

- #4220 IUI / #4220 IUI Memo

5 boxes of 25 individually sterile packed, disposable devices.

- #4220 Plus

3 boxes of 25 individually sterile packed, disposable devices.

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