**#4219 EMTRAC A**

Embryo Transfer Catheter

**HOW SUPPLIED:** #4219 Emtrac A

Important product information, read carefully before use.

**Disposable – discard after single use****Re-use may cause sexually transmitted infections****Caution: federal (USA) law restricts this device to sale by or on the order of a physician****Sterile: unless pouch is damaged or opened****DESCRIPTION:**

The device is an embryo transfer catheter. It is a one piece catheter with markings to indicate depth, has a handle with luer lock and a positioner for depth adjustment. Intended patient population: mature women

**INDICATIONS FOR USE:**

Introduction of embryos into the uterine cavity.

**CONTRAINDICATIONS:**

Do not use in the presence of cervical infection, after recent pelvic inflammatory disease, sexually transmitted disease or who is currently pregnant. Do not use if the package is opened/damaged; discard immediately.

**WARNINGS:**

Please refer to published medical literature before performing this procedure. The device should be used by healthcare professionals only. Forcing the catheter may result in damage to the endometrial tissue and bleeding. If passage of the device through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**PRECAUTIONS:**

By virtue of its flexibility and small diameter, the device can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroflexion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

**ADVERSE REACTIONS:**

No adverse reactions associated with the use of the device have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping and superficial haemorrhage may be expected. In every case the directions for use should be followed; taking note of contra-indications, warnings and precautions.

1639

**Manufactured by:**

Gynetics Medical Products N.V.  
Rembert Dodoensstraat 51  
3920 Lommel, Belgium

**PACKAGED:**

5 boxes of 25 individually sterile packed, disposable devices.

**DIRECTIONS FOR USE:**

1. Check the expiry date of the device
2. Position the patient in the lithotomy position. Insert vaginal speculum to expose cervix. Determine the depth and direction of the patient's uterus.
3. Set the positioner as previously measured.
4. Attach the syringe to the luer lock and load the embryo(s) into the device, following your standard protocol.
5. Gently advance the catheter into the external cervical os and advance it through the cervical canal into the uterine cavity until resistance is felt. The distal tip of the catheter should then be at the uterine fundus. Gently rotate the catheter to remove any excess mucus from the tip.

**Note:** If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. The device should not be used in case of resistance to achieve insertion.

6. Slowly press the plunger of the syringe and expel the embryos into the uterine cavity.
7. Following the complete expulsion of the embryos, slowly withdraw the catheter from the uterus. Inspect the catheter to verify embryo release
8. Remove the speculum and allow the patient to rest.
9. Discard the catheter after use in accordance with local regulation.

Catalogue number

Batch code

Manufacturer

Temperature limit

Consult instructions for use

CE marking of conformity

Use-by date

Sterilized using irradiation

Do not use if package is damaged

Do not reuse

Medical device

Do not resterilize

Single sterile barrier system

Single sterile barrier system with protective packaging

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