#4219 EMTRAC C-200 / #4219 EMTRAC C-240

Embryo Transfer catheter

HOW SUPPLIED:**REF** #4219 Emtrac C-200**REF** #4219 Emtrac C-240

Important information, read carefully before use

Disposable – discard after single use**Re-use may cause sexually transmitted infections****Caution: federal (USA) law restricts this device to sale by or on the order of a physician****Sterile: unless pouch is damaged or opened****DESCRIPTION:**

The device is an embryo transfer catheter. It is a single use loading catheter and has a handle with luer lock. The device is available in 2 sizes. To be used in combination with a re-usable guiding catheter. Intended patient population: mature women undergoing medical assisted reproduction procedures

INDICATIONS FOR USE:

Introduction of embryos into the uterine cavity.

CONTRAINDICATIONS:

Do not use in the presence of cervical infection, after recent pelvic inflammatory disease, sexually transmitted disease or who is currently pregnant. Do not use if the package is opened/damaged; discard immediately.

**WARNINGS:**

Please refer to published medical literature before performing this procedure. The device should be used by healthcare professionals only trained in medical assisted reproduction procedures

Forcing the catheter may result in damage to the endometrial tissue and bleeding. If passage of the device through the cervical canal is not found to be easily achievable, the device should not be forced and the possibility of the presence of pathologic cervical stenosis considered. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS:

By virtue of its flexibility and small diameter, the device can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken, prior to its insertion, to ascertain the depth of the uterus and any uterine anteversion or retroversion present in individual patients. (This is preferably done during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized).

ADVERSE REACTIONS:

No adverse reactions associated with the use of the device have been reported. However, as with any instrument, which passes through the internal cervical os, mild cramping and superficial haemorrhage may be expected. In every case the direction for use should be followed; taking note of contra-indications, warnings and precautions.

1639

**Manufactured by:**

Gynetics Medical Products N.V.,
Rembert Dodoensstraat 51
3920 Lommel, Belgium

PACKAGED:

3 boxes of 25 individually sterile packed, disposable devices.

DIRECTIONS FOR USE:

The device should be used in combination with a re-usable guiding catheter.

1. Check the expiry date of the device.
2. Position the patient in the lithotomy position. Insert vaginal speculum to expose cervix. Determine the depth and direction of the patient's uterus.
3. Gently advance the guiding catheter into the external cervical os and advance it through the cervical canal into the uterine cavity until resistance is felt. Gently rotate the catheter to remove any excess mucus from the tip.

Note: If any difficulty with insertion or bending of the catheter occurs, gently rotate the catheter while advancing it close to the cervix. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. The device should not be used in case of resistance to achieve insertion.

4. Attach the syringe to the luer lock of the loading catheter. Load the embryo(s) into the loading catheter following your standard protocol. Insert the loading catheter in the guiding catheter.
5. Slowly press the plunger of the syringe and expel the embryos into the uterine cavity.
6. Following the complete expulsion of the embryos, slowly withdraw the catheters from the uterus. Inspect the loading catheter to verify embryo release.
7. Remove the speculum and allow the patient to rest.
8. Discard the catheter in accordance with local regulation.

Catalogue number

Batch code

Manufacturer

Temperature limit

Consult instructions for use

CE marking of conformity

Use-by date

Sterilized using irradiation

Do not use if package if damaged

Do not reuse

Medical device

Do not resterilize

Single sterile barrier system

Single sterile barrier system with protective packaging

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