

SEMTRAC SET 4.2

SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

FEATURES

- soft bullet shaped tip to minimize trauma
- three piece catheter for difficult transfers
- 5 free mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- reinforced shaft of loading catheter for optimal control
- handle with luer lock for a more secure fitting to syringes
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

LOADING CATHETER TIP OPTIONS



CLASSIC

TIP-IN-SIGHT

GYNÓ INFO SEMTRAC SET 4.2

Guiding catheter (A) length	158 mm
tube OD	2.55 mm
Loading catheter (C) length	232 mm
tube OD	1.70 mm
total volume	0.23 ml
total volume minus hub	0.081 ml

Sterile ETO sterilized


Packaging 75/shipping box

Quality control MEA and LAL

Order numbers

- Semtrac Set 4.2 #2000 Set 4.2
- Semtrac Set 4.2 Tip-in-Sight #2000 Set 4.2 TIS

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.



GMP Belgium Semtrac Set 4.2 06/2018



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be



Embryo Transfer