

SEMTRAC SET 5

SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

FEATURES

- soft bullet shaped tip to minimize trauma
- three piece catheter for difficult transfers
- 5 **free** mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- reinforced shaft of loading catheter for optimal control
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter
- **NEW:** available with **Tip-in-Sight (TIS)** technology for precise positioning of the catheter tip under ultrasound guidance

LOADING CATHETER TIP OPTIONS



CLASSIC

TIP-IN-SIGHT

GYNO INFO SEMTRAC SET 5

Guiding catheter (A) length 178 mm
tube OD 2.55 mm

Loading catheter (C) length 252 mm
tube OD 1.70 mm
total volume 0.22 ml
total volume minus hub 0.089 ml

Sterile ETO sterilized


Packaging 75/shipping box

Quality control MEA and LAL

Order numbers

- Semtrac Set 5 #2000 Set 5 ☒
- Semtrac Set 5 Trial #2000 Set 5 T ☒
- Semtrac Set 5 Tip-in-Sight #2000 Set 5 TIS ☒

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.



Guiding catheter (A)

Loading catheter (C) #2000 Set 5

Loading catheter (C) #2000 Set 5 TIS

Mandrel (B)

GMP Belgium Semtrac Set 5 06/2018

