

SEMTRAC SET W-180

SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

FEATURES

- soft bullet shaped tip to minimize trauma
- 5 free mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

LOADING CATHETER TIP OPTIONS



GYNÓ INFO SEMTRAC SET W-180

Guiding catheter (A) length 141 mm
tube OD 2.55 mm

Loading catheter (C) length 214 mm
tube OD 1.70 mm
total volume 0.23 ml
total volume minus hub 0.074 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

Order numbers

- Semtrac Set W-180 #2000 W-180 ☒
- Semtrac Set W-180 Trial #2000 W-180 T ☒
- Semtrac Set W-180 Tip-in-Sight #2000 W-180 TIS ☒

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.

