

IUI

UTERINE INSEMINATOR FOR WASHED SPERMATOZOA

The IUI is an intrauterine insemination cannula with two lateral distal ports for optimal sperm distribution.

FEATURES


- positioner to prevent damage to fundus
- two lateral ports at distal end for optimal sperm distribution
- arrows on handle to indicate lateral ports at distal end
- markings on the outer tube
- handle with luer lock for a more secure fitting to syringes



ORIGINAL PICTURE

GYNO INFO IUI

Length	206 mm (± 2 mm)
Tube OD	1.85 mm (± 0.10 mm)
Sterile	sterilized by gamma irradiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	#4220 



GMP Belgium IUI 09/2014



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine