

IUI PLUS

UTERINE INSEMINATOR FOR WASHED SPERMATOZOA WITH MANDREL

The IUI Plus is a set for difficult cervical passage, consisting of two parts:

- A) Loading cannula
- B) Mandrel

FEATURES


- positioner for depth adjustment
- a combination of a soft, flexible intrauterine cannula and a solid cervix cannula
- markings on outer tube for easy positioning
- mandrel to prevent mucus to enter the loading cannula
- handle with luer lock for a more secure fitting to syringes



ORIGINAL PICTURE

GYNO INFO IUI PLUS

Length	A) 216 mm (\pm 3 mm) B) 240 mm (\pm 2 mm)
Tube OD	A) 1.60 mm (\pm 0.10 mm) B) 1.00 mm (\pm 0.05 mm)
Total volume	A) 0.42 ml
Total volume minus hub	A) 0.258 ml
Sterile	sterilized by gamma irradiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4220 Plus



GMP Belgium IUI Plus 09/2014



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine