

# SMOOZE IUI

## FLEXIBLE INTRAUTERINE INSEMINATION CANNULA

The Smooze IUI is a flexible intrauterine insemination cannula with a smooth tip for less traumatic inseminations.

### FEATURES

- smooth tip with one distal side hole
- very flexible and therefore less traumatic
- handle with luer lock for a more secure fitting to syringes
- available in 3 lengths:

- A) Short      7 cm (not including handle)  
B) Medium    14 cm (not including handle)  
C) Long       20 cm (not including handle)



ORIGINAL PICTURE

### GYNO INFO SMOOZE IUI

Length :            A) 100 mm ( $\pm$  2 mm)  
                          B) 169 mm ( $\pm$  2 mm)  
                          C) 230 mm ( $\pm$  2 mm)




Tube OD :           2 mm ( $\pm$  0.05 mm)


Volume :            A) min. 0.030 ml  
                          B) min. 0.057 ml  
                          C) min. 0.080 ml

Sterile :              ETO sterilized

Packaging :         A) 125/shipping box  
                          B) 125/shipping box  
                          C) 75/shipping box

Quality control :   MEA and LAL

Order number :    A) #4225   
                          B) #4225 M   
                          C) #4225 L 



GMP Belgium Smooze IUI 09/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine