

# EMTRAC A

## EMBRYO TRANSFER CATHETER WITH POSITIONER

The Emtrac A is a combination of a soft, flexible intrauterine catheter and a solid cervix catheter.


### FEATURES

- positioner for depth adjustment
- a combination of a soft, flexible intrauterine catheter and a solid cervix catheter
- markings on outer tube for easy positioning
- handle with luer lock for a more secure fitting to syringes
- also implemented for intrauterine insemination, when a thinner cannula is required

ORIGINAL PICTURE

### GYNO INFO EMTRAC A

Length	216 mm (± 3 mm)
Tube OD	1.60 mm (± 0.10 mm)
Total volume	0.42 ml
Total volume minus hub	0.258 ml
Sterile	sterilized by radiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	<b>#4219 A</b> ⓘ



GMP Belgium Emtrac A 04/2011



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine