

SUPERIOR FAS

FOLLICLE ASPIRATION SET FEATURES & ACCESSORIES

NEEDLE TYPES

3 needle types, available in various sizes and dimensions to determine your own specification:

Single Lumen Single Lumen Flush Double Lumen Flush

NEEDLE TIPS

2 types of needle tips to guarantee clean facet cutting:

Back Cut Bevel  Lancet Bevel 



TIP MARKINGS

2 types of tip markings to offer excellent ultrasound visibility of the needle tip:

Dual Laser Marking  Embossed Marking 

HANDLES

2 types of handles available for individual comfort (only for single lumen flush & double lumen flush):

U shaped handle  Y shaped handle 

TRANSLUCENT TUBING


3 types of translucent tubing in any desired length (maximum length thick wall tubing is 120cm/47.2 inch):

Regular Tubing for standard oocyte aspiration set up
Soft Tubing allows control of the fluid
Thick Wall Tubing no kinking of the tubing



CORKS


2 types of silicone corks for a secure connection to test tubes:

Small 14x18x20 mm 0.55 x 0.71 x 0.79 inch
Medium 21x 27x30 mm 0.83 x 1.06 x 1.18 inch



OTHER ACCESSORIES - SPECIALS

3-way Valve  Luer Lock Needle (For direct connection to a syringe) 

Male/male Connector  Stainless steel flushing cannula 



Contact us for other or special enquiries!