

# GYNOSAMPLER

## ENDOMETRIAL HISTOLOGY

The Gynosampler is the American name for the Probet which is sold in the USA and Canada. The Gynosampler is an endometrial suction curette for a histologic lining biopsy.

### FEATURES

- double headed piston secures a steady vacuum
- tissue collection through side hole
- smooth rounded distal tip and small size for minimal patient discomfort
- depth markings for ease in locating depth of insertion in the uterus



ORIGINAL PICTURE

### GYNO INFO GYNOSAMPLER

Length	262 mm (± 3 mm)
Tube OD	2.93 mm (± 0.03 mm)
Sterile	ETO sterilized
Packaging	500/shipping box (Canada only) 250/shipping box (US only)
Order number	<b>#9164-CAN</b> ☺ (Canada only) <b>#9164</b> ☺ (US only)



GMP Belgium Gynosampler 04/2017



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine