

HANDLING PIPETTE

PIPETTE FOR HANDLING EMBRYOS AND OOCYTES

The handling pipette is a single-use IVF laboratory micropipette used for transferring oocytes and embryos from one dish to another.

FEATURES

- IVF laboratory micropipette for dish-to-dish transfer of oocytes and embryos
- will not scratch dishes
- available in three lengths: 45, 55 and 94 mm
- available in two packaging variations: single packed and packed by ten pieces



ORIGINAL PICTURE

GYNO INFO HANDLING PIPETTE

Length 45, 55 or 94 mm (\pm 3 mm)


Tube OD 1.25 mm (\pm 0.05 mm)
Tube ID 1.00 mm (\pm 0.05 mm)

Sterile sterilized by gamma irradiation

Packaging 250/shipping box (1/pouch)
1250/shipping box (10/pouch)

Quality control MEA and LAL

Order number #3145 (45 mm, 1/pouch)
#3100 (55 mm, 1/pouch)
#3125 (94 mm, 1/pouch)
#3145 TP (45 mm, 10/pouch)
#3100 TP (55 mm, 10/pouch)
#3125 TP (94 mm, 10/pouch)



GMP Belgium Handling Pipette 09/2016

