

STRIP-ID

PIPETTE FOR OOCYTES MANIPULATION

The Strip-ID is a single-use IVF laboratory micropipette used for oocytes manipulation and handling of embryos and oocytes

FEATURES

- superior quality that ensures smooth edges for precise denudation
- strict individual microscopic quality inspection
- flexible, unbreakable material that will not scratch dishes
- convenient colour-coded multi-packaging that allows easy release of a single pipette
- extensive variety of sizes available to suit your application
- fits the "Pipet-ID(s)" and is compatible with most popular used micropipettors



ORIGINAL PICTURE

ORIG. PIC. 200X
ENLARGED

GYNO INFO STRIP-ID

Length	90 mm (\pm 3 mm)		
Inner diameter μ m	125 μ m	145 μ m	200 μ m
	130 μ m	150 μ m	275 μ m
	135 μ m	170 μ m	300 μ m
	140 μ m	175 μ m	
Sterile	sterilized by gamma irradiation		
Packaging	10 pipettes / vial 15 vials / inner box 75 vials of 10 pipettes / shipping box		
Quality control	MEA and LAL		
Order number	#3000 Strip-ID125/10  #3000 Strip-ID130/10  #3000 Strip-ID135/10  #3000 Strip-ID140/10  #3000 Strip-ID145/10  #3000 Strip-ID150/10  #3000 Strip-ID170/10  #3000 Strip-ID175/10  #3000 Strip-ID200/10  #3000 Strip-ID275/10  #3000 Strip-ID300/10 		

GMP Belgium Strip-ID 12/2016



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine

GYNO INFO

STRIP-ID

Order number	ID	Colour
#3000 Strip-ID125/10	125 µm	blue
#3000 Strip-ID130/10	130 µm	yellow
#3000 Strip-ID135/10	135 µm	yellow
#3000 Strip-ID140/10	140 µm	white
#3000 Strip-ID145/10	145 µm	white
#3000 Strip-ID150/10	150 µm	black
#3000 Strip-ID170/10	170 µm	red
#3000 Strip-ID175/10	175 µm	red
#3000 Strip-ID200/10	200 µm	dark red
#3000 Strip-ID275/10	275 µm	green
#3000 Strip-ID300/10	300 µm	transparant



Easy release of a single pipette



Colour-coded multi-packaging

