

# PIPET-ID

## MICROPIPETTOR

The Pipet-ID and Pipet-IDs are reusable IVF laboratory micropipettors developed to fit the "Strip-ID" micropipettes and most other popular micropipettes (Ø 0.9 mm OD)

### FEATURES

- Available in two different sizes to fit every size of hand
  - **Pipet-ID** for regular to large hands
  - **Pipet-IDs** for small to regular hands
- Ergonomic design
- Adjustable aspiration volume (min ≈ 0.2 µl, max ≈ 3 µl)
- Total expulsion feature
- Stainless steel table holder included
- Spare parts kit included (includes extra seals, plunger wire and hexagon key)
- Designed for use with the "Strip-ID" micropipettes. Also fits most plastic micropipettes (Ø 0.9 mm OD, 0.6 mm ID)



ORIGINAL PICTURE

### GYNO INFO PIPET-ID

Length (setting on position 3) 143 mm

Weight 57 gr


Order number #3000 Pipet-ID

### GYNO INFO PIPET-IDs

Length (setting on position 3) 120 mm

Weight 43 gr

Order number #3000 Pipet-IDs



GMP Belgium Pipet-ID 11/2016



Gynétics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNÉTICS**<sup>®</sup>  
Practise better Medicine