

# SEMTRAC C

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac C is a soft embryo transfer catheter.

### FEATURES

- soft bullet shaped tip to minimize trauma
- one piece economical catheter
- markings to indicate depth
- metal tube inside for better control
- handle with luer lock for a more secure fitting to syringes
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

### TIP OPTIONS



### GYNO INFO SEMTRAC C

Length	252 mm	
Tube OD	1.70 mm	
Total volume	0.20 ml	
Total volume minus hub	0.065 ml	
Sterile	ETO sterilized	
Packaging	75/shipping box	
Quality control	MEA and LAL	
Order numbers		
- Semtrac C	#2000 C	ⓧ
- Semtrac C Tip-in-Sight	#2000 C TIS	ⓧ

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.

