

SEMTRAC SET W-240

SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

FEATURES

- soft bullet shaped tip to minimize trauma
- 5 free mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

LOADING CATHETER TIP OPTIONS

CLASSIC

TIP-IN-SIGHT

GYNÓ INFO SEMTRAC SET W-240

Guiding catheter (A) length	203 mm
tube OD	2.55 mm
Loading catheter (C) length	277 mm
tube OD	1.70 mm
total volume	0.26 ml
total volume minus hub	0.098 ml

Sterile ETO sterilized


Packaging 75/shipping box

Quality control MEA and LAL

Order numbers

- Semtrac Set W-240 #2000 W-240 ②
- Semtrac Set W-240 Trial #2000 W-240 T ②
- Semtrac Set W-240 Tip-in-Sight #2000 W-240 TIS ②

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.



GMP Belgium Semtrac W-240 06/2018

