

THE TULIP MEMO

EMBRYO TRANSFER CATHETER WITH TULIP SHAPED BULB

The Tulip is a set for difficult cervical passage, consisting of two parts: a guiding catheter (A) and a loading catheter (C).

FEATURES

- the Tulip tip enables easy atraumatic entry through the cervical canal
- mandrel available for extra difficult transfers
- positioner for depth adjustment
- a combination of a soft, flexible intrauterine catheter and a solid cervix catheter with memory
- markings on outer tube for easy positioning
- memo function of the guiding catheter
- 5 free mandrels to preshape the guiding catheter in every box of 25 pieces
- handle with luer lock fitting for a more secure connection to syringes
- reinforced metal support of loading catheter for optimal control



ORIGINAL PICTURE

GYNO INFO THE TULIP MEMO

Length A) 200 mm
C) 266 mm

Tube OD A) 2.20 mm
C) 1.00 mm

Total volume C) 0.17 ml
Total volume minus hub C) 0.035 ml


Tulip bulb 3 mm

Sterile sterilized by radiation

Packaging 75/shipping box

Quality control MEA and LAL

Order number #4000 Memo ②



GMP Belgium The Tulip Memo 06/2018