

EMTRAC DELPHIN

EMBRYO TRANSFER CATHETER WITH POSITIONER AND SOFT UTERINE CATHETER

The combination of a soft flexible intrauterine catheter and a solid cervix catheter creates an universal catheter for easy insertion in the most difficult cervical tract.

FEATURES

- positioner to prevent damage to fundus
- one piece catheter
- markings to indicate depth
- extra soft uterine catheter to ensure safe insertion through cervix
- handle with luer lock for a more secure fitting to syringes
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

TIP OPTIONS



GYNO INFO EMTRAC DELPHIN

Length	210 mm
Tube OD	2.00 mm
Total volume	0.23 ml
Total volume minus hub	0.072 ml
Sterile	ETO sterilized
Packaging	75/shipping box
Quality control	MEA and LAL

Order numbers

- Emtrac Delphin #4219 Delphin ②
- Emtrac Delphin Tip-in-Sight #4219 Delphin TIS ②

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.

