


RENEWED EMBRYO TRANSFER CATHETER PORTFOLIO

The renewed Embryo Transfer catheters offer you a variety of possibilities to enhance your practice



The Tulip

EMBRYO TRANSFER CATHETER SETS WITH BULB SHAPED TIP

The renewed Tulip Embryo Transfer catheter range offers a variety of possibilities to support you enhancing your practice.

- The **bulb shaped** guiding catheters assist atraumatic entry into the cervical canal and are available in different shapes.
- The **Memo function** allows shaping the guiding catheter in any desired silhouette to suit various patients' anatomies.
- The Tulip **loading catheters** are available with a Classic tip or a Tip-in-Sight (TIS), which supports precise positioning under ultrasound guidance.

Guiding Catheter bulb shape options



The Tulip S
Oval shaped bulb



The Tulip L
Rigid round shaped bulb



The Tulip LS
SOFT round shaped bulb

Loading Catheter tip options





Tip-in-Sight (TIS) Echogenic



CHOOSE YOUR TULIP SET

| PRODUCT | | ORDER NUMBER | FEATURES | | | | |
|-------------------------|--------------------------------|--------------------|---------------|--------------------------------|--------------------|--|--|
| | | | Echogenic tip | Preshaped with curve indicator | Memory function | Mandrel available 5 per 25 sets | |
| | THE TULIP S | #4000 Tulip S | | ✓ | | | |
| | THE TULIP S Memo | #4000 S Memo | | | ✓ | ✓ | |
| | THE TULIP S Tip-in-Sight | #4000 Tulip S TIS | ✓ | ✓ | | | |
| Oval shaped bulb | THE TULIP S Memo Tip-in-Sight | #4000 S Memo TIS | ✓ | | ✓ | ✓ | |
| | THE TULIP L | #4000 Tulip L | | ✓ | | | |
| | THE TULIP L Memo | #4000 L Memo | | | ✓ | ✓ | |
| | THE TULIP L Tip-in-Sight | #4000 Tulip L TIS | ✓ | ✓ | | | |
| Rigid round shaped bulb | THE TULIP L Memo Tip-in-Sight | #4000 L Memo TIS | ✓ | | ✓ | ✓ | |
| | THE TULIP LS | #4000 Tulip LS | | ✓ | | | |
| | THE TULIP LS Memo | #4000 LS Memo | | | ✓ | ✓ | |
| | THE TULIP LS Tip-in-Sight | #4000 Tulip LS TIS | ✓ | ✓ | | | |
| SOFT round shaped bulb | THE TULIP LS Memo Tip-in-Sight | #4000 LS Memo TIS | ✓ | | ✓ | ✓ | |

GYNO INFO

| Length | guiding catheter: 200 mm loading catheter: 266 mm | | | |
|--|---|--|--|--|
| Tube OD | guiding catheter: 2.20 mm loading catheter: 1.00 mm | | | |
| Total volume Total volume minus hub | loading catheter: 0.21 ml loading catheter: 0.069 ml | | | |
| Bulb shape | oval shaped: 3 mm rigid/soft round shaped: 2.3 mm | | | |
| Sterile | sterilized by gamma irradiation | | | |
| Packaging | 75/shipping box | | | |
| Quality control | MEA and LAL | | | |



Emtrac Delphin / Semtrac / The Tulip

EMBRYO TRANSFER CATHETERS WITH TIP-IN-SIGHT (TIS) TECHNOLOGY

The TIS technology supports precise positioning of the loading catheter tip under ultrasound guidance. The new range of Embryo Transfer catheters with the TIS technology is an extension of the existing portfolio with Classic loading tips.



GYNO INFO

| The Tulip TIS EMBRYO TRANSFER CATHETER SETS WITH BULB SHAPED TIP | | | | | | | | | | |
|--|---------------------|------------------|---------|------------------|---------|---------|--------------------|---------------|--|--|
| PRODUCT | ORDER NUMBER | GUIDING CATHETER | | LOADING CATHETER | | VOLUME | | | | |
| | | LENGTH | TUBE OD | LENGTH | TUBE OD | TOTAL | TOTAL MINUS HUB | BULB SHAPE | | |
| THE TULIP S Tip-in-Sight | #4000 Tulip S TIS | 200 mm | 2.20 mm | 266 mm | 1.00 mm | 0.21 ml | 0.069 ml | 3 mm | | |
| THE TULIP S Memo Tip-in-Sight | #4000 S Memo TIS | 200 mm | 2.20 mm | 266 mm | 1.00 mm | 0.21 ml | 0.069 ml | 3 mm | | |
| THE TULIP L Tip-in-Sight | #4000 Tulip L TIS | 200 mm | 2.20 mm | 266 mm | 1.00 mm | 0.21 ml | 0.069 ml | 2.3 mm | | |
| THE TULIP L Memo Tip-in-Sight | #4000 L Memo TIS | 200 mm | 2.20 mm | 266 mm | 1.00 mm | 0.21 ml | 0.069 ml | 2.3 mm | | |
| THE TULIP LS Tip-in-Sight | #4000 Tulip LS TIS | 200 mm | 2.20 mm | 266 mm | 1.00 mm | 0.21 ml | 0.069 ml | 2.3 mm | | |
| THE TULIP LS Memo Tip-in-Sight | #4000 LS Memo TIS | 200 mm | 2.20 mm | 266 mm | 1.00 mm | 0.21 ml | 0.069 ml | 2.3 mm | | |
| Semtrac TIS SOFT EMBRYO TRANSFER CATHETER | | | | | | | | | | |
| PRODUCT | ORDER NUMBER | GUIDING CATHETER | | LOADING CATHETER | | VOLUME | | | | |
| | | LENGTH | TUBE OD | LENGTH | TUBE OD | TOTAL | TOTAL MINUS HUB | | | |
| Semtrac C Tip-in-Sight | #2000 C TIS | - | - | 252 mm | 1.70 mm | 0.20 ml | 0.065 ml | | | |
| Semtrac Set 4.2 Tip-in-Sight | #2000 Set 4.2 TIS | 158 mm | 2.55 mm | 232 mm | 1.70 mm | 0.23 ml | 0.081 ml | | | |
| Semtrac Set 5 Tip-in-Sight | #2000 Set 5 TIS | 178 mm | 2.55 mm | 252 mm | 1.70 mm | 0.22 ml | 0.089 ml | | | |
| Semtrac Set W-180 Tip-in-Sight | #2000 Set W-180 TIS | 141 mm | 2.55 mm | 214 mm | 1.70 mm | 0.23 ml | 0.074 ml | | | |
| Semtrac Set W-240 Tip-in-Sight | #2000 Set W-240 TIS | 203 mm | 2.55 mm | 277 mm | 1.70 mm | 0.26 ml | 0.098 ml | | | |
| Emtrac Delphin TIS EMBRYO TRANSFER CATHETER WITH POSITIONER | | | | | | | | | | |
| PRODUCT | ORDER NUMBER | GUIDING CATHETER | | LOADING CATHETER | | VOLUME | | | | |
| | | LENGTH | TUBE OD | LENGTH | TUBE OD | TOTAL | TOTAL MINUS HUB | | | |
| Emtrac Delphin Tip-in-Sight | #4219 Delphin TIS | - | - | 210 mm | 2.00 mm | 0.23 ml | 0.072 ml | | | |

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.

90240118 001



