

# GYNO INFO

## INTRAUTERINE INSEMINATION CANNULA RANGE

PRODUCT	ORDER NUMBER
IUI	#4220
IUI Memo	#4220 Memo
IUI Plus	#4220 Plus
Smooze IUI	#4225
Smooze IUI Medium	#4225 M
Smooze IUI Long	#4225 L
Seminor Rigid	#4502 A
Seminor Flexible	#4502 B

GMP Belgium IUI Range 09/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# IUI

## UTERINE INSEMINATOR FOR WASHED SPERMATOZOA

The IUI is an intrauterine insemination cannula with two lateral distal ports for optimal sperm distribution.

### FEATURES


- positioner to prevent damage to fundus
- two lateral ports at distal end for optimal sperm distribution
- arrows on handle to indicate lateral ports at distal end
- markings on the outer tube
- handle with luer lock for a more secure fitting to syringes



ORIGINAL PICTURE

### GYNO INFO IUI

Length	206 mm (± 2 mm)
Tube OD	1.85 mm (± 0.10 mm)
Sterile	sterilized by gamma irradiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	#4220 



GMP Belgium IUI 09/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# IUI MEMO

## UTERINE INSEMINATOR FOR WASHED SPERMATOZOA WITH MEMORY FUNCTION

The IUI Memo is an intrauterine insemination cannula with MEMO function, which proactively enables the distal tip to bend according to cervical flexion.


### FEATURES


- positioner to prevent damage to fundus
- two lateral ports at distal end for optimal sperm distribution
- arrows on handle to indicate lateral ports at distal end
- markings on outer tube
- memory function
- handle with luer lock for a more secure fitting to syringes



ORIGINAL PICTURE

### GYNO INFO IUI MEMO

Length	206 mm (± 2 mm)
Tube OD	1.85 mm (± 0.10 mm)
Sterile	sterilized by gamma irradiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	#4220 Memo 



GMP Belgium IUI Memo 09/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# IUI PLUS

## UTERINE INSEMINATOR FOR WASHED SPERMATOZOA WITH MANDREL

The IUI Plus is a set for difficult cervical passage, consisting of two parts:

- A) Loading cannula
- B) Mandrel


### FEATURES

- positioner for depth adjustment
- a combination of a soft, flexible intrauterine cannula and a solid cervix cannula
- markings on outer tube for easy positioning
- mandrel to prevent mucus to enter the loading cannula
- handle with luer lock for a more secure fitting to syringes

ORIGINAL PICTURE

### GYNO INFO IUI PLUS

Length	A) 216 mm ( $\pm$ 3 mm) B) 240 mm ( $\pm$ 2 mm)
Tube OD	A) 1.60 mm ( $\pm$ 0.10 mm) B) 1.00 mm ( $\pm$ 0.05 mm)
Total volume	A) 0.42 ml
Total volume minus hub	A) 0.258 ml
Sterile	sterilized by gamma irradiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4220 Plus 



GMP Belgium IUI Plus 09/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# SMOOZE IUI

## FLEXIBLE INTRAUTERINE INSEMINATION CANNULA

The Smooze IUI is a flexible intrauterine insemination cannula with a smooth tip for less traumatic inseminations.

### FEATURES

- smooth tip with one distal side hole
- very flexible and therefore less traumatic
- handle with luer lock for a more secure fitting to syringes
- available in 3 lengths:

- A) Short 7 cm (not including handle)  
B) Medium 14 cm (not including handle)  
C) Long 20 cm (not including handle)



ORIGINAL PICTURE

### GYNO INFO SMOOZE IUI

Length : A) 100 mm ( $\pm$  2 mm)  
B) 169 mm ( $\pm$  2 mm)  
C) 230 mm ( $\pm$  2 mm)




Tube OD : 2 mm ( $\pm$  0.05 mm)


Volume : A) min. 0.030 ml  
B) min. 0.057 ml  
C) min. 0.080 ml

Sterile : ETO sterilized

Packaging : A) 125/shipping box  
B) 125/shipping box  
C) 75/shipping box

Quality control : MEA and LAL

Order number : A) #4225   
B) #4225 M   
C) #4225 L 



GMP Belgium Smooze IUI 09/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# SEMINOR RIGID/FLEXIBLE

## MULTI-PURPOSE CANNULA FOR ASSISTED HUMAN REPRODUCTION

- A) The Seminor Rigid is a semi-rigid cannula for intrauterine and cervical insemination.
- B) The Seminor Flexible is a cannula for intrauterine and cervical insemination.


### FEATURES

- the Seminor Rigid (A) can be shaped according to the degree of uterine flexion
- the Seminor Flexible (B) enables easy insertion with minimal patient discomfort
- the conical proximal end closes the cervix adequately for maximum sperm transfer
- as result of a direct connection between Seminor and syringe, the loss of spermatozoa during the transfer is minimized




### GYNO INFO SEMINOR RIGID/FLEXIBLE

Length	48 mm ( $\pm$ 0.5 mm)
Volume	min. 0.098 ml
Sterile	sterilized by gamma irradiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	A) #4502 A  B) #4502 B 



GMP Belgium Seminor 09/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# GYNO INFO

## EMBRYO TRANSFER CATHETER RANGE

### EMBRYO TRANSFER CATHETERS FOR EASY TRANSFER

PRODUCT	ORDER NUMBER
Emtrac A	#4219 A
Emtrac Delphin	#4219 Delphin
Emtrac Delphin Tip-in-Sight	#4219 Delphin TIS
Semtrac C	#2000 C
Semtrac C Tip-in-Sight	#2000 C TIS

### EMBRYO TRANSFER CATHETER SETS FOR DIFFICULT TRANSFER

PRODUCT	ORDER NUMBER
Emtrac Set	#4219 Set
Semtrac Set 4.2	#2000 Set 4.2
Semtrac Set 4.2 Tip-in-Sight	#2000 Set 4.2 TIS
Semtrac Set 5	#2000 Set 5
Semtrac Set 5 Trial	#2000 Set 5 T
Semtrac Set 4.2 Tip-in-Sight	#2000 Set 4.2 TIS
Semtrac Set W-180	#2000 W-180
Semtrac Set W-180 Trial	#2000 W-180 T
Semtrac Set W-180 Tip-in-Sight	#2000 W-180 TIS
Semtrac Set W-240	#2000 W-240
Semtrac Set W-240 Trial	#2000 W-240 T
Semtrac Set W-180 Tip-in-Sight	#2000 W-180 TIS

GMP Belgium ET Range 05/2014



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# GYNO INFO

## EMBRYO TRANSFER CATHETER RANGE

### EMBRYO TRANSFER CATHETER SETS FOR DIFFICULT TRANSFER

PRODUCT	ORDER NUMBER
The Tulip	#4000 The Tulip
The Tulip Memo	#4000 Memo
The Tulip S	#4000 Tulip S
The Tulip S Memo	#4000 Tulip S Memo
The Tulip S Tip-in-Sight	#4000 Tulip S TIS
The Tulip S Memo Tip-in-Sight	#4000 Tulip S Memo TIS
The Tulip L	#4000 Tulip L
The Tulip L Memo	#4000 Tulip L Memo
The Tulip L Tip-in-Sight	#4000 Tulip L TIS
The Tulip L Memo Tip-in-Sight	#4000 Tulip L Memo TIS
The Tulip LS	#4000 Tulip LS
The Tulip LS Memo	#4000 Tulip LS Memo
The Tulip LS Tip-in-Sight	#4000 Tulip LS TIS
The Tulip LS Memo Tip-in-Sight	#4000 Tulip LS Memo TIS

### REUSABLE EMBRYO TRANSFER CATHETERS FOR DIFFICULT TRANSFER

PRODUCT	ORDER NUMBER
Guiding Catheter 200	#Gyreflex-200FC
Guiding Catheter 240	#Gyreflex-240FC

GMP Belgium ET Range 05/2014





# EMTRAC A

## EMBRYO TRANSFER CATHETER WITH POSITIONER

The Emtrac A is a combination of a soft, flexible intrauterine catheter and a solid cervix catheter.


### FEATURES

- positioner for depth adjustment
- a combination of a soft, flexible intrauterine catheter and a solid cervix catheter
- markings on outer tube for easy positioning
- handle with luer lock for a more secure fitting to syringes
- also implemented for intrauterine insemination, when a thinner cannula is required

ORIGINAL PICTURE

### GYNO INFO EMTRAC A

Length	216 mm (± 3 mm)
Tube OD	1.60 mm (± 0.10 mm)
Total volume	0.42 ml
Total volume minus hub	0.258 ml
Sterile	sterilized by radiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	#4219 A 



GMP Belgium Emtrac A 04/2011



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# EMTRAC DELPHIN

## EMBRYO TRANSFER CATHETER WITH POSITIONER AND SOFT UTERINE CATHETER

The combination of a soft flexible intrauterine catheter and a solid cervix catheter creates an universal catheter for easy insertion in the most difficult cervical tract.

### FEATURES

- positioner to prevent damage to fundus
- one piece catheter
- markings to indicate depth
- extra soft uterine catheter to ensure safe insertion through cervix
- handle with luer lock for a more secure fitting to syringes
- **NEW:** available with **Tip-in-Sight (TIS) technology** for precise positioning of the catheter tip under ultrasound guidance

### TIP OPTIONS



### GYNO INFO EMTRAC DELPHIN

Length	210 mm
Tube OD	2.00 mm
Total volume	0.23 ml
Total volume minus hub	0.072 ml
Sterile	ETO sterilized
Packaging	75/shipping box
Quality control	MEA and LAL

#### Order numbers

- Emtrac Delphin #4219 Delphin ②
- Emtrac Delphin Tip-in-Sight #4219 Delphin TIS ②

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.



# SEMTRAC C

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac C is a soft embryo transfer catheter.

### FEATURES

- soft bullet shaped tip to minimize trauma
- one piece economical catheter
- markings to indicate depth
- metal tube inside for better control
- handle with luer lock for a more secure fitting to syringes
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

### TIP OPTIONS



### GYNO INFO SEMTRAC C

Length	252 mm	
Tube OD	1.70 mm	
Total volume	0.20 ml	
Total volume minus hub	0.065 ml	
Sterile	ETO sterilized	
Packaging	75/shipping box	
Quality control	MEA and LAL	
Order numbers		
- Semtrac C	#2000 C	ⓧ
- Semtrac C Tip-in-Sight	#2000 C TIS	ⓧ

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.



# EMTRAC SET

## EMBRYO TRANSFER CATHETER WITH POSITIONER

The Emtrac Set is a set for difficult cervical passage, consisting of three parts:


- A) Guiding catheter
- B) Mandrel
- C) Loading catheter


### FEATURES

- positioner for depth adjustment
- a combination of a soft, flexible intrauterine catheter and a solid cervix catheter
- markings on outer tube for easy positioning
- soft loading catheter to minimize trauma
- handle with luer lock for a more secure fitting to syringes

ORIGINAL PICTURE

### GYNO INFO EMTRAC SET

Length	A) 216 mm B) 240 mm C) 266 mm
Tube OD	A) 1.60 mm B) 1.00 mm C) 1.00 mm
Total volume	C) 0.20 ml
Total volume minus hub	C) 0.069 ml
Sterile	sterilized by radiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4219 Set 



GMP Belgium Emtrac Set 04/2011



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# SEMTRAC SET 4.2

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

### FEATURES

- soft bullet shaped tip to minimize trauma
- three piece catheter for difficult transfers
- 5 free mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- reinforced shaft of loading catheter for optimal control
- handle with luer lock for a more secure fitting to syringes
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

### LOADING CATHETER TIP OPTIONS



CLASSIC

TIP-IN-SIGHT

### GYNO INFO SEMTRAC SET 4.2

Guiding catheter (A) length	158 mm
tube OD	2.55 mm
Loading catheter (C) length	232 mm
tube OD	1.70 mm
total volume	0.23 ml
total volume minus hub	0.081 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

#### Order numbers

- Semtrac Set 4.2 #2000 Set 4.2
- Semtrac Set 4.2 Tip-in-Sight #2000 Set 4.2 TIS

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.



# SEMTRAC SET 5

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

### FEATURES

- soft bullet shaped tip to minimize trauma
- three piece catheter for difficult transfers
- 5 free mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- reinforced shaft of loading catheter for optimal control
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

### LOADING CATHETER TIP OPTIONS



### GYNO INFO SEMTRAC SET 5

Guiding catheter (A) length 178 mm  
tube OD 2.55 mm

Loading catheter (C) length 252 mm  
tube OD 1.70 mm  
total volume 0.22 ml  
total volume minus hub 0.089 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

#### Order numbers

- Semtrac Set 5 #2000 Set 5 ☒
- Semtrac Set 5 Trial #2000 Set 5 T ☒
- Semtrac Set 5 Tip-in-Sight #2000 Set 5 TIS ☒

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.



# SEMTRAC SET W-180

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

### FEATURES

- soft bullet shaped tip to minimize trauma
- 5 free mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

### LOADING CATHETER TIP OPTIONS



### GYNÓ INFO SEMTRAC SET W-180

Guiding catheter (A) length 141 mm  
tube OD 2.55 mm

Loading catheter (C) length 214 mm  
tube OD 1.70 mm  
total volume 0.23 ml  
total volume minus hub 0.074 ml

Sterile ETO sterilized


Packaging 75/shipping box

Quality control MEA and LAL

#### Order numbers

- Semtrac Set W-180 #2000 W-180 ☒
- Semtrac Set W-180 Trial #2000 W-180 T ☒
- Semtrac Set W-180 Tip-in-Sight #2000 W-180 TIS ☒

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.





# SEMTRAC SET W-240

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

### FEATURES

- soft bullet shaped tip to minimize trauma
- 5 free mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter
- **NEW:** available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance

### LOADING CATHETER TIP OPTIONS

CLASSIC

TIP-IN-SIGHT

### GYNÓ INFO SEMTRAC SET W-240

Guiding catheter (A) length	203 mm
tube OD	2.55 mm
Loading catheter (C) length	277 mm
tube OD	1.70 mm
total volume	0.26 ml
total volume minus hub	0.098 ml

Sterile ETO sterilized


Packaging 75/shipping box

Quality control MEA and LAL

#### Order numbers

- Semtrac Set W-240 #2000 W-240 ②
- Semtrac Set W-240 Trial #2000 W-240 T ②
- Semtrac Set W-240 Tip-in-Sight #2000 W-240 TIS ②

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.






# THE TULIP

## PRESHAPED EMBRYO TRANSFER CATHETER WITH TULIP SHAPED BULB

The Tulip is a set for difficult cervical passage, consisting of two parts: a guiding catheter (A) and a loading catheter (C).

### FEATURES


- preshaped guiding catheter, ideally curved for easy entry when uterus is flexed
- the Tulip tip enables easy atraumatic entry through the cervical canal
- positioner for depth adjustment
- markings on outer tube for easy positioning
- handle with luer lock fitting for a more secure connection to syringes
- reinforced metal support of loading catheter for optimal control
- arrow marking on outer tube for curve indication



ORIGINAL PICTURE

### GYNO INFO THE TULIP

Length	A) 200 mm C) 266 mm
Tube OD	A) 2.20 mm C) 1.00 mm
Total volume	C) 0.17 ml
Total volume minus hub	C) 0.035 ml
Tulip bulb	3 mm
Sterile	sterilized by radiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4000 The Tulip <sup>®</sup>



GMP Belgium The Tulip 04/2011



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine


# THE TULIP MEMO

## EMBRYO TRANSFER CATHETER WITH TULIP SHAPED BULB

The Tulip is a set for difficult cervical passage, consisting of two parts: a guiding catheter (A) and a loading catheter (C).

### FEATURES


- the Tulip tip enables easy atraumatic entry through the cervical canal
- mandrel available for extra difficult transfers
- positioner for depth adjustment
- a combination of a soft, flexible intrauterine catheter and a solid cervix catheter with memory
- markings on outer tube for easy positioning
- memo function of the guiding catheter
- 5 free mandrels to preshape the guiding catheter in every box of 25 pieces
- handle with luer lock fitting for a more secure connection to syringes
- reinforced metal support of loading catheter for optimal control



ORIGINAL PICTURE

### GYNO INFO THE TULIP MEMO

Length	A) 200 mm C) 266 mm
Tube OD	A) 2.20 mm C) 1.00 mm
Total volume	C) 0.17 ml
Total volume minus hub	C) 0.035 ml
Tulip bulb	3 mm
Sterile	sterilized by radiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4000 Memo ②




GMP Belgium The Tulip Memo 04/2011

## EMBRYO TRANSFER CATHETER SETS WITH BULB SHAPED TIP

The renewed Tulip Embryo Transfer catheter range offers a variety of possibilities to support you enhancing your practice.

- The **bulb shaped** guiding catheters assist atraumatic entry into the cervical canal and are available in different shapes.
- The **Memo function** allows shaping the guiding catheter in any desired silhouette to suit various patients' anatomies.
- The Tulip **loading catheters** are available with a Classic tip or a Tip-in-Sight (TIS), which supports precise positioning under ultrasound guidance.

### Guiding Catheter **bulb shape** options



**The Tulip S**  
Oval shaped bulb




**The Tulip L**  
Rigid round shaped bulb



**The Tulip LS**  
**SOFT** round shaped bulb

### Loading Catheter **tip** options



**Classic**



**Tip-in-Sight (TIS)**  
Echogenic




## CHOOSE YOUR TULIP SET

PRODUCT	ORDER NUMBER	FEATURES			
		Echogenic tip	Preshaped with curve indicator	Memory function	Mandrel available 5 per 25 sets
 Oval shaped bulb	THE TULIP S	#4000 Tulip S		✓	
	THE TULIP S Memo	#4000 S Memo			✓
	THE TULIP S Tip-in-Sight	#4000 Tulip S TIS	✓	✓	
	THE TULIP S Memo Tip-in-Sight	#4000 S Memo TIS	✓		✓
 Rigid round shaped bulb	THE TULIP L	#4000 Tulip L		✓	
	THE TULIP L Memo	#4000 L Memo			✓
	THE TULIP L Tip-in-Sight	#4000 Tulip L TIS	✓	✓	
	THE TULIP L Memo Tip-in-Sight	#4000 L Memo TIS	✓		✓
 <b>SOFT</b> round shaped bulb	THE TULIP LS	#4000 Tulip LS		✓	
	THE TULIP LS Memo	#4000 LS Memo			✓
	THE TULIP LS Tip-in-Sight	#4000 Tulip LS TIS	✓	✓	
	THE TULIP LS Memo Tip-in-Sight	#4000 LS Memo TIS	✓		✓

## GYNO INFO

Length	guiding catheter: 200 mm loading catheter: 266 mm
Tube OD	guiding catheter: 2.20 mm loading catheter: 1.00 mm
Total volume Total volume minus hub	loading catheter: 0.21 ml loading catheter: 0.069 ml
Bulb shape	oval shaped: 3 mm rigid/soft round shaped: 2.3 mm
Sterile	sterilized by gamma irradiation
Packaging	75/shipping box
Quality control	MEA and LAL



# GUIDING CATHETER

## EMBRYO TRANSFER CATHETER FOR DIFFICULT EMBRYO TRANSFERS

The Guiding Catheter is a reusable, economical and soft embryo transfer catheter for difficult cervical passage.

### FEATURES

- visible markings for exact positioning
- combined with a single use loading catheter: Emtrac C-200 or -240
- bulb tip enables easy atraumatic cervical passage
- flexible end which enables it to follow the cervical canal
- echogenic for optimal visualisation during ultrasound guided transfers




ORIGINAL PICTURE

### GYNO INFO GUIDING CATHETER 200 & 240

Length	200 - 240 mm
Tube ID	1.10 mm
Tube OD	2.50/1.90 mm
Packaging	1/shipping box
Order number	#Gyneflex-200FC #Gyneflex-240FC

### GYNO INFO EMTRAC C-200 & C-240

Length	286 - 297 mm (± 2 mm)
Tube ID	0.60 mm (± 0.05 mm)
OD	1.00 mm (± 0.05 mm)
Total volume	0.207 ml - 0.210 ml
Total volume minus hub	0.075 ml - 0.078 ml
Sterile	sterilized by gamma irradiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4219 C-200 ☺ #4219 C-240 ☺



GMP Belgium Guiding Catheter 05/2014




Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# GYNETICS®

Practise better Medicine

## FOLLICLE ASPIRATION NEEDLES





Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# GYNO INFO

## SUPERIOR FOLLICLE ASPIRATION NEEDLE RANGE

PRODUCT	PAGE
Features and accessories	4
Single lumen	5-6
Single lumen flush	7-8
Double lumen flush	9-10



## PREMIUM FOLLICLE ASPIRATION NEEDLE RANGE

PRODUCT	PAGE
Single lumen	11-12





# SUPERIOR FAS

## FOLLICLE ASPIRATION SET FEATURES & ACCESSORIES



### NEEDLE TYPES

3 needle types, available in various sizes and dimensions to determine your own specification:

Single Lumen      Single Lumen Flush      Double Lumen Flush



### NEEDLE TIPS

2 types of needle tips to guarantee clean facet cutting:

Back Cut Bevel       Lancet Bevel 

### TIP MARKINGS

2 types of tip markings to offer excellent ultrasound visibility of the needle tip:

Dual Laser Marking       Embossed Marking 

### HANDLES

2 types of handles available for individual comfort (only for single lumen flush & double lumen flush):

U shaped handle       Y shaped handle 

### TRANSLUCENT TUBING


3 types of translucent tubing in any desired length (maximum length thick wall tubing is 120cm/47.2 inch):

Regular Tubing      for standard oocyte aspiration set up  
Soft Tubing      allows control of the fluid  
Thick Wall Tubing      no kinking of the tubing



### CORKS

2 types of silicone corks for a secure connection to test tubes:

Small      14x18x20 mm      0.55 x 0.71 x 0.79 inch  
Medium      21x 27x30 mm      0.83 x 1.06 x 1.18 inch



### OTHER ACCESSORIES - SPECIALS

3-way Valve       Luer Lock Needle (For direct connection to a syringe) 

Male/male Connector       Stainless steel flushing cannula 



Contact us for other or special enquiries!

# SUPERIOR FAS SINGLE LUMEN

## CUSTOM-MADE FOLLICLE ASPIRATION SET

The Superior FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.


### FEATURES


- sharp lancet and back cut bevel tip available that ensure a clean facet cutting
- dual laser and embossed echogenic marked tip available which offer stable and excellent vision during ultrasound
- various needle sizes for ovaries with complicated and/or limited access
- translucent tubing in the same size as needle for safe oocyte retrieval
- arrow marking on handle for bevel point indication
- handle shape facilitates rotation
- custom-made



ORIGINAL PICTURE


### GYNO INFO SUPERIOR FAS SINGLE LUMEN

Sterile	sterilized by gamma irradiation
Packaging	50/shipping box (standard configuration)
Quality control	MEA and LAL
Order number	<b>#4551</b>  (complete the configurator and receive your order number)



# SUPERIOR FAS SINGLE LUMEN

## CONFIGURE YOUR OWN CUSTOM-MADE SUPERIOR FOLLICLE ASPIRATION SET IN 6 STEPS



COPY OR DOWNLOAD THIS CONFIGURATOR AT [WWW.GYNETICS.BE](http://WWW.GYNETICS.BE),  
COMPLETE AND RETURN FOR INQUIRY TO RECEIVE YOUR CUSTOMIZED ORDER NUMBER

### STEP 1 CHOOSE YOUR NEEDLE

TICK BOX

Gauge	OD/ID (mm)	OD/ID (inch)	Length (mm)	Length (inch)	Control	
16G	1.6/1.1	0.063/0.043	400	15.75	Super rigid	<input type="checkbox"/>
16G	1.6/1.3	0.063/0.051	400	15.75	Medium	<input type="checkbox"/>
16G	1.6/1.3	0.063/0.051	340	13.39	Medium	<input type="checkbox"/>
17G	1.5/0.9	0.059/0.035	400	15.75	Super rigid	<input type="checkbox"/>
17G	1.5/0.9	0.059/0.035	350	13.78	Super rigid	<input type="checkbox"/>
17G	1.4/1.0	0.055/0.039	400	15.75	Rigid	<input type="checkbox"/>
17G	1.4/1.0	0.055/0.039	350	13.78	Rigid	<input type="checkbox"/>
17G	1.4/1.0	0.055/0.039	330	12.99	Rigid	<input type="checkbox"/>
17G	1.4/1.2	0.055/0.047	400	15.75	Flexible	<input type="checkbox"/>
17G	1.4/1.2	0.055/0.047	340	13.39	Flexible	<input type="checkbox"/>
18G	1.2/1.0	0.047/0.039	400	15.75	Medium	<input type="checkbox"/>
18G	1.2/1.0	0.047/0.039	350	13.78	Medium	<input type="checkbox"/>
19G	1.0/0.65	0.039/0.026	330	12.99	Medium	<input type="checkbox"/>
20G	0.9/0.6	0.035/0.024	370	14.57	Medium	<input type="checkbox"/>

The working length of each needle is the length minus 50 mm/1.968 inch

### STEP 2 CHOOSE YOUR TUBING FROM NEEDLE TO CORK

Type	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>
No tubing, Female Luer Connector at end of needle (Luer Lock Needle)		step 3 to 5 n/a	<input type="checkbox"/>

### STEP 3 CHOOSE YOUR SILICONE CORK

Type	Size (mm)	Size (inch)	
Small Cork	14 x 18 x 20	0.55 x 0.71 x 0.79	<input type="checkbox"/>
Medium Cork	21 x 27 x 30	0.83 x 1.06 x 1.18	<input type="checkbox"/>





### STEP 4 CHOOSE YOUR TUBING FROM CORK TO PUMP

Type	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>

### STEP 5 CHOOSE YOUR CONNECTION TO THE PUMP

Female Luer Connector (standard)	<input type="checkbox"/>
Male Luer Connector	<input type="checkbox"/>
No Connector	<input type="checkbox"/>

### STEP 6 CHOOSE YOUR VARIATION OF THE TIP

Back cut bevel with embossed echogenic marking		<input type="checkbox"/>
Lancet bevel with embossed echogenic marking		<input type="checkbox"/>
Back cut bevel with dual laser echogenic marking		<input type="checkbox"/>
Lancet bevel with dual laser echogenic marking		<input type="checkbox"/>

# SUPERIOR FAS SINGLE LUMEN FLUSH

## CUSTOM-MADE FOLLICLE ASPIRATION SET

The Superior FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.


### FEATURES


- flushing connector enables flushing and aspirating during complicated oocyte retrieval or when few mature oocytes occur
- flushing and aspirating procedure is performed separately through one lumen
- sharp lancet and back cut bevel tip available that ensure a clean facet cutting
- dual laser and embossed echogenic marked tip available which offer stable and excellent vision during ultrasound
- translucent tubing in the same size as needle for safe oocyte retrieval
- various needle sizes for ovaries with complicated and/or limited access
- handle shape facilitates rotation
- custom-made



ORIGINAL PICTURE

### GYNO INFO SUPERIOR FAS SINGLE LUMEN FLUSH


Sterile	sterilized by gamma irradiation
Packaging	50/shipping box (standard configuration)
Quality control	MEA and LAL
Order number	#4551  (complete the configurator and receive your order number)



# SUPERIOR FAS SINGLE LUMEN FLUSH

## CONFIGURE YOUR OWN CUSTOM-MADE SUPERIOR FOLLICLE ASPIRATION SET IN 7 STEPS

COPY OR DOWNLOAD THIS CONFIGURATOR AT [WWW.GYNETICS.BE](http://WWW.GYNETICS.BE),  
COMPLETE AND RETURN FOR INQUIRY TO RECEIVE YOUR CUSTOMIZED ORDER NUMBER





### STEP 1 CHOOSE YOUR NEEDLE AND HANDLE

TICK BOX

Gauge	OD/ID (mm)	OD/ID (inch)	Length (mm)	Length (inch)	Control	
17G	1.40/1.20	0.055/0.047	352	13.858	Flexible	<input type="checkbox"/>
17G	1.40/1.00	0.055/0.039	352	13.858	Rigid	<input type="checkbox"/>

The working length of each needle is 300 mm/11.811 inch



### STEP 2 CHOOSE YOUR FLUSH TUBING

Tubing	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>

### STEP 3 CHOOSE YOUR TUBING FROM NEEDLE TO CORK

Tubing	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>

### STEP 4 CHOOSE YOUR SILICONE CORK

Type	Size (mm)	Size (inch)	
Small Cork	14 x 18 x 20	0.55 x 0.71 x 0.79	<input type="checkbox"/>
Medium Cork	21 x 27 x 30	0.83 x 1.06 x 1.18	<input type="checkbox"/>



### STEP 5 CHOOSE YOUR TUBING FROM CORK TO PUMP

Type	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>

### STEP 6 CHOOSE YOUR CONNECTION TO THE PUMP

Female Luer Connector (standard)	<input type="checkbox"/>
Male Luer Connector	<input type="checkbox"/>
No Connector	<input type="checkbox"/>

### STEP 7 CHOOSE YOUR VARIATION OF THE TIP

Back cut bevel with embossed echogenic marking		<input type="checkbox"/>
Lancet bevel with embossed echogenic marking		<input type="checkbox"/>
Back cut bevel with dual laser echogenic marking		<input type="checkbox"/>
Lancet bevel with dual laser echogenic marking		<input type="checkbox"/>



# SUPERIOR FAS DOUBLE LUMEN FLUSH

## CUSTOM-MADE FOLLICLE ASPIRATION SET

The Superior FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.

### FEATURES


- flushing connector enables flushing and aspirating during complicated oocyte retrieval or when few mature oocytes occur
- double lumen allows simultaneous flushing and aspirating
- dual laser and embossed echogenic marked tip available which offer stable and excellent vision during ultrasound
- translucent tubing in the same size as needle for safe oocyte retrieval
- various needle sizes for ovaries with complicated and/or limited access
- handle shape facilitates rotation
- custom-made



ORIGINAL PICTURE

### GYNO INFO SUPERIOR FAS DOUBLE LUMEN FLUSH


Sterile	sterilized by gamma irradiation
Packaging	50/shipping box (standard configuration)
Quality control	MEA and LAL
Order number	<b>#4551</b> ② (complete the configurator and receive your order number)



# SUPERIOR FAS DOUBLE LUMEN FLUSH

CONFIGURE YOUR OWN  
CUSTOM-MADE SUPERIOR  
FOLLICLE ASPIRATION SET IN 7 STEPS

COPY OR DOWNLOAD THIS CONFIGURATOR AT [WWW.GYNETICS.BE](http://WWW.GYNETICS.BE),  
COMPLETE AND RETURN FOR INQUIRY TO RECEIVE YOUR CUSTOMIZED ORDER NUMBER





**STEP 1 CHOOSE YOUR NEEDLE AND HANDLE**

TICK BOX

Gauge	OD/ID (mm)	OD/ID (inch)	Length (mm)	Length (inch)	Control	
16G	outer lumen 1.65/1.40 inner lumen 1.20/1.00	0.065/0.055 0.047/0.039	352	13.858	Rigid	<input type="checkbox"/>
17G	outer lumen 1.50/1.20 inner lumen 1.00/0.80	0.059/0.047 0.039/0.032	352	13.858	Rigid	<input type="checkbox"/>

The working length of each needle is 300 mm/11.811 inch



**STEP 2 CHOOSE YOUR FLUSH TUBING**

Tubing	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>

**STEP 3 CHOOSE YOUR TUBING FROM NEEDLE TO CORK**

Tubing	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>

**STEP 4 CHOOSE YOUR SILICONE CORK**

Type	Size (mm)	Size (inch)	
Small Cork	14 x 18 x 20	0.55 x 0.71 x 0.79	<input type="checkbox"/>
Medium Cork	21 x 27 x 30	0.83 x 1.06 x 1.18	<input type="checkbox"/>


**STEP 5 CHOOSE YOUR TUBING FROM CORK TO PUMP**

Type	Length needed	Remarks	
Regular Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Soft Tubing	_____ cm/inch	n/a	<input type="checkbox"/>
Thick Wall Tubing	_____ cm/inch	max. 120 cm/47.2 inch	<input type="checkbox"/>

**STEP 6 CHOOSE YOUR CONNECTION TO THE PUMP**

Female Luer Connector (standard)	<input type="checkbox"/>
Male Luer Connector	<input type="checkbox"/>
No Connector	<input type="checkbox"/>

**STEP 7 CHOOSE YOUR VARIATION OF THE TIP**

Dual laser echogenic marking 	<input type="checkbox"/>
Embossed echogenic marking 	<input type="checkbox"/>

# PREMIUM FAS SINGLE LUMEN

## FOLLICLE ASPIRATION SET

The Premium FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.

### FEATURES

- standard follicle aspiration needle with regular tubing
- sharp lancet and back cut bevel tip available that ensure a clean facet cutting
- diamond laser ultrasound markings on tip which offer stable and excellent vision during ultrasound
- translucent tubing in the same size as needle for safe oocyte retrieval
- arrow marking on handle for bevel point indication
- handle shape facilitates rotation



ORIGINAL PICTURE


### GYNO INFO PREMIUM FAS SINGLE LUMEN

Sterile                   sterilized by gamma irradiation

Packaging             50/shipping box (standard configuration)

Quality control       MEA and LAL

- |              |                 |   |
|--------------|-----------------|---|
| Order number | #4551 NS-AS1    | Ⓜ |
|              | #4551 NS-BS1    | Ⓜ |
|              | #4551 NS-CS1    | Ⓜ |
|              | #4551 NS-DS1    | Ⓜ |
|              | #4551 NS-AS1BCB | Ⓜ |
|              | #4551 NS-CS1BCB | Ⓜ |
|              | #4551 NS-DS1BCB | Ⓜ |






# GYNO INFO

## PREMIUM FAS SINGLE LUMEN

① Order number	Gauge	OD/ID (mm)	Length (mm)	Working length (mm)	Bevel	Control
#4551 NS-AS1	17G	1.4/1.2	400	350	Lancet	Flexible
#4551 NS-BS1	17G	1.4/1.2	340	290	Lancet	Flexible
#4551 NS-CS1	17G	1.4/1.0	400	350	Lancet	Rigid
#4551 NS-DS1	17G	1.4/1.0	340	290	Lancet	Rigid
#4551 NS-AS1BCB	17G	1.4/1.2	400	350	Back cut bevel	Flexible
#4551 NS-CS1BCB	17G	1.4/1.0	400	350	Back cut bevel	Rigid
#4551 NS-DS1BCB	17G	1.4/1.0	340	290	Back cut bevel	Rigid

Order number	Gauge	OD/ID (inch)	Length (inch)	Working length (inch)	Bevel	Control
#4551 NS-AS1	17G	0.055/0.047	15.75	13.779	Lancet	Flexible
#4551 NS-BS1	17G	0.055/0.047	13.39	11.417	Lancet	Flexible
#4551 NS-CS1	17G	0.055/0.039	15.75	13.779	Lancet	Rigid
#4551 NS-DS1	17G	0.055/0.039	13.39	11.417	Lancet	Rigid
#4551 NS-AS1BCB	17G	0.055/0.047	15.75	13.779	Back cut bevel	Flexible
#4551 NS-CS1BCB	17G	0.055/0.039	15.75	13.779	Back cut bevel	Rigid
#4551 NS-DS1BCB	17G	0.055/0.039	13.39	11.417	Back cut bevel	Rigid

- ② Tubing from needle to cork : Regular tubing 100 cm / 39.369 inch
- ③ Size silicone cork : Small 14 x 18 x 20 mm / 0.551 x 0.709 x 0.787 inch
- ④ Tubing from cork to pump : Regular tubing 100 cm / 39.369 inch



# GYNO INFO

## SAMPLING RANGE

PRODUCT ORDER NUMBER

Detecta #4163

Probet #4164

Gynosampler #9164

In/out #4165

Echosampler #5164

GMP Belgium Sampling 11/2017



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# DETECTA

## CERVICAL SCREENING DEVICE DESIGNED TO OBTAIN OPTIMAL SAMPLES

The Detecta is a simple but efficient cervical sampling device.

### FEATURES

- allows in-office sampling
- systematic brushing of the transitional zone
- helps to detect neoplasia at an early stage
- easy to handle



ORIGINAL PICTURE

### GYNO INFO DETECTA

Length	194 mm
Packaging	1250/shipping box
Order number	#4163 ☒

GMP Belgium Detecta 06/09





Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# MANUAL DETECTA

- gently insert the Detecta into the cervix
- rotate only once to obtain the sample
- gently withdraw the Detecta without rotating
- smear sample on a slide with one rotation of the Detecta
- fixate the slide immediately
- sample is ready for examination




# PROBET

## ENDOMETRIAL HISTOLOGY

The Probet is an endometrial suction curette for a histologic lining biopsy.

### FEATURES

- double headed piston secures a steady vacuum
- tissue collection through side hole
- smooth rounded distal tip and small size for minimal patient discomfort
- depth markings for ease in locating depth of insertion in the uterus



ORIGINAL PICTURE

### GYNO INFO PROBET

Length	262 mm (± 3 mm)
Tube OD	2.93 mm (± 0.03 mm)
Sterile	ETO sterilized
Packaging	75/shipping box
Order number	<b>#4164</b> ⓘ



GMP Belgium Probet 04/2011




Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be


**GYNETICS**<sup>®</sup>  
Practise better Medicine

# MANUAL PROBET


1) Insert Probet through cervical canal while piston is fully engaged in outer tube ("sheath") until wall contact is felt, taking into account hystero-metry and curvature.




2) Hold sheath. Pull piston back completely until stopper is reached, in order to create maximum negative pressure within sheath.



3) Roll sheath between fingers 4 times – while simultaneously moving it back and forth inside the uterus – for comprehensive sampling.



4) Remove device. Cut distal tip just below side opening. Push piston back into sheath to expell sample in transport medium.




# GYNOSAMPLER

## ENDOMETRIAL HISTOLOGY

The Gynosampler is the American name for the Probet which is sold in the USA and Canada. The Gynosampler is an endometrial suction curette for a histologic lining biopsy.

### FEATURES


- double headed piston secures a steady vacuum
- tissue collection through side hole
- smooth rounded distal tip and small size for minimal patient discomfort
- depth markings for ease in locating depth of insertion in the uterus



ORIGINAL PICTURE

### GYN INFO GYNOSAMPLER

Length	262 mm (± 3 mm)
Tube OD	2.93 mm (± 0.03 mm)
Sterile	ETO sterilized
Packaging	500/shipping box (Canada only) 250/shipping box (US only)
Order number	<b>#9164-CAN</b> ☺ (Canada only) <b>#9164</b> ☺ (US only)



GMP Belgium Gynosampler 04/2017



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# IN/OUT

## ENDOMETRIAL HISTOLOGY

The IN/OUT is an endometrial suction curette for a histologic lining biopsy.

### FEATURES

- double headed piston secures a steady vacuum
- tissue collection through the tip
- smooth rounded distal tip and small size for minimal patient discomfort
- depth markings for ease in locating depth of insertion in the uterus



ORIGINAL PICTURE

### GYNO INFO IN/OUT

Length	262 mm (± 3 mm)
Tube OD	2.93 mm (± 0.03 mm)
Sterile	ETO sterilized
Packaging	75/shipping box
Order number	<b>#4165</b> 



GMP Belgium IN/OUT 04/2011




Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine




# MANUAL IN / OUT


1) Insert In/out through cervical canal while piston is fully engaged in outer tube ('sheath') until wall contact is felt, taking into account hysteroemetry and curvature.




2) Hold sheath. Pull piston back completely until stopper is reached, in order to create maximum negative pressure within sheath.



3) Roll sheath between fingers 4 times – while simultaneously moving it back and forth inside the uterus – for comprehensive sampling.



4) Remove device. Push piston back into sheath to expell sample in transport medium.



# ECHOSAMPLER

## INTRACAVITAL PATHOLOGY AND ENDOMETRIAL HISTOLOGY

The EchoSampler is a flexible sampling catheter used for dual diagnosis: intra cavital pathology and endometrial histology.

### FEATURES

- dual diagnosis: cavital pathology and endometrial histology
- for ultrasound guided biopsy
- for saline infusion sonohysterography
- implemented in case of abnormal uterine bleeding
- flexible sampling catheter
- depth markings for ease in locating depth of insertion in the uterus
- tissue collection and saline infusion through side hole
- smooth rounded distal tip and small size for minimal patient discomfort and easy passage through cervical canal
- handle with luer lock for a more secure fitting to syringes
- echogenic line to facilitate visualisation when using ultrasound



ORIGINAL PICTURE

### GYNO INFO ECHOSAMPLER

Length	180 mm (± 3 mm)
Tube OD	2.93 mm (± 0.04 mm)
Sterile	ETO sterilized
Packaging	75/shipping box
Order number	#5164 ☒



GMP Belgium Echosampler 02/2017



# GYNO INFO

## IVF LABORATORY RANGE

### PRODUCT

### ORDER NUMBER

Handling Pipette

#3145  
#3100  
#3125  
#3145 TP  
#3100 TP  
#3125 TP

Strip-ID

#3000 Strip-ID125/10  
#3000 Strip-ID130/10  
#3000 Strip-ID135/10  
#3000 Strip-ID140/10  
#3000 Strip-ID145/10  
#3000 Strip-ID150/10  
#3000 Strip-ID170/10  
#3000 Strip-ID175/10  
#3000 Strip-ID200/10  
#3000 Strip-ID275/10  
#3000 Strip-ID300/10

Pipet-ID

#3000 Pipet-ID  
#3000 Pipet-IDs

GMP Belgium IVF Laboratory 09/2016



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# HANDLING PIPETTE

## PIPETTE FOR HANDLING EMBRYOS AND OOCYTES

The handling pipette is a single-use IVF laboratory micropipette used for transferring oocytes and embryos from one dish to another.

### FEATURES

- IVF laboratory micropipette for dish-to-dish transfer of oocytes and embryos
- will not scratch dishes
- available in three lengths: 45, 55 and 94 mm
- available in two packaging variations: single packed and packed by ten pieces



ORIGINAL PICTURE

### GYNO INFO HANDLING PIPETTE







Length 45, 55 or 94 mm ( $\pm$  3 mm)


Tube OD 1.25 mm ( $\pm$  0.05 mm)  
Tube ID 1.00 mm ( $\pm$  0.05 mm)

Sterile sterilized by gamma irradiation

Packaging 250/shipping box (1/pouch)  
1250/shipping box (10/pouch)

Quality control MEA and LAL

Order number **#3145**  (45 mm, 1/pouch)  
**#3100**  (55 mm, 1/pouch)  
**#3125**  (94 mm, 1/pouch)  
**#3145 TP**  (45 mm, 10/pouch)  
**#3100 TP**  (55 mm, 10/pouch)  
**#3125 TP**  (94 mm, 10/pouch)



GMP Belgium Handling Pipette 09/2016


# PIPET-ID

## MICROPIPETTOR

The Pipet-ID and Pipet-IDs are reusable IVF laboratory micropipettors developed to fit the "Strip-ID" micropipettes and most other popular micropipettes (Ø 0.9 mm OD)

### FEATURES

- Available in two different sizes to fit every size of hand
  - **Pipet-ID** for regular to large hands
  - **Pipet-IDs** for small to regular hands
- Ergonomic design
- Adjustable aspiration volume (min ≈ 0.2 µl, max ≈ 3 µl)
- Total expulsion feature
- Stainless steel table holder included
- Spare parts kit included (includes extra seals, plunger wire and hexagon key)
- Designed for use with the "Strip-ID" micropipettes. Also fits most plastic micropipettes (Ø 0.9 mm OD, 0.6 mm ID)



ORIGINAL PICTURE

### GYNO INFO PIPET-ID

Length (setting on position 3) 143 mm

Weight 57 gr


Order number #3000 Pipet-ID

### GYNO INFO PIPET-IDs

Length (setting on position 3) 120 mm

Weight 43 gr

Order number #3000 Pipet-IDs



GMP Belgium Pipet-ID 10/2017



Colour differences may occur between batches due to variations in the colouring process and be considered as normal. GMP will not accept any claim or remark concerning these eventual differences.



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNÉTIQS**<sup>®</sup>  
Practise better Medicine

# STRIP-ID

## PIPETTE FOR OOCYTES MANIPULATION

The Strip-ID is a single-use IVF laboratory micropipette used for oocytes manipulation and handling of embryos and oocytes

### FEATURES

- superior quality that ensures smooth edges for precise denudation
- strict individual microscopic quality inspection
- flexible, unbreakable material that will not scratch dishes
- convenient colour-coded multi-packaging that allows easy release of a single pipette
- extensive variety of sizes available to suit your application
- fits the "Pipet-ID(s)" and is compatible with most popular used micropipettors

			
ORIGINAL PICTURE		ORIG. PIC. 200X ENLARGED	
<b>GYN INFO STRIP-ID</b>			
Length	90 mm (± 3 mm)		
Inner diameter µm	125 µm	145 µm	200 µm
	130 µm	150 µm	275 µm
	135 µm	170 µm	300 µm
	140 µm	175 µm	
Sterile	sterilized by gamma irradiation		
Packaging	10 pipettes / vial 15 vials / inner box 75 vials of 10 pipettes / shipping box		
Quality control	MEA and LAL		
Order number	#3000 Strip-ID125/10  #3000 Strip-ID130/10  #3000 Strip-ID135/10  #3000 Strip-ID140/10  #3000 Strip-ID145/10  #3000 Strip-ID150/10  #3000 Strip-ID170/10  #3000 Strip-ID175/10  #3000 Strip-ID200/10  #3000 Strip-ID275/10  #3000 Strip-ID300/10 		

GMP Belgium Strip-ID 01/2018



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine

# GYNO INFO

## STRIP-ID

Order number	ID	Colour
#3000 Strip-ID125/10	125 µm	blue
#3000 Strip-ID130/10	130 µm	yellow
#3000 Strip-ID135/10	135 µm	yellow
#3000 Strip-ID140/10	140 µm	white
#3000 Strip-ID145/10	145 µm	white
#3000 Strip-ID150/10	150 µm	black
#3000 Strip-ID170/10	170 µm	red
#3000 Strip-ID175/10	175 µm	red
#3000 Strip-ID200/10	200 µm	dark red
#3000 Strip-ID275/10	275 µm	green
#3000 Strip-ID300/10	300 µm	transparant



Easy release of a single pipette



Colour-coded multi-packaging

