


SUPERIOR FAS DOUBLE LUMEN FLUSH

CUSTOM-MADE FOLLICLE ASPIRATION SET

The Superior FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.

FEATURES

- flushing connector enables flushing and aspirating during complicated oocyte retrieval or when few mature oocytes occur
- double lumen allows simultaneous flushing and aspirating
- dual laser and embossed echogenic marked tip available which offer stable and excellent vision during ultrasound
- translucent tubing in the same size as needle for safe oocyte retrieval
- various needle sizes for ovaries with complicated and/or limited access
- handle shape facilitates rotation
- custom-made



ORIGINAL PICTURE

GYNO INFO SUPERIOR FAS DOUBLE LUMEN FLUSH

Sterile	sterilized by gamma irradiation
Packaging	50/shipping box (standard configuration)
Quality control	MEA and LAL
Order number	#4551 ② (complete the configurator and receive your order number)



Gynetics Medical Products n.v.
 Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
 tel. +32 11 64 58 72 / fax +32 11 64 06 77
 info@gynetics.be / www.gynetics.be

GYNETICS[®]
 Practise better Medicine