

# SUPERIOR FAS

## FOLLICLE ASPIRATION SET

### NEEDLE TYPES

3 needle types, available in various sizes and dimensions to determine your own specification:

- Single Lumen
- Single Lumen Flush
- Double Lumen Flush

### NEEDLE TIPS

2 types of needle tips to guarantee clean facet cutting:

- Back Cut Bevel 
- Lancet Bevel 

### TIP MARKINGS

2 types of tip markings to offer excellent ultrasound visibility of the needle tip:

- Dual Laser Marking 
- Embossed Marking 

### HANDLES

2 types of handles available for individual comfort (only for single lumen flush & double lumen flush):

- U shaped handle 
- Y shaped handle 



### TRANSLUCENT TUBING

3 types of translucent tubing in any desired length (maximum length thick wall tubing is 120cm/47.2 inch):

- Regular Tubing for standard oocyte aspiration set up
- Soft Tubing allows control of the fluid
- Thick Wall Tubing no kinking of the tubing

### CORKS

2 types of silicone corks for a secure connection to test tubes:

Small	14x18x20 mm	0.55 x 0.71 x 0.79 inch		
Medium	21x27x30 mm	0.83 x 1.06 x 1.18 inch		



Contact our sales department for other or special enquiries!