

SUPERIOR FAS

FOLLICLE ASPIRATION SET

NEEDLE TYPES

3 needle types, available in various sizes and dimensions to determine your own specification:

Single Lumen Single Lumen Flush Double Lumen Flush

NEEDLE TIPS

2 types of needle tips to guarantee clean facet cutting:

Back Cut Bevel  Lancet Bevel 

TIP MARKINGS

2 types of tip markings to offer excellent ultrasound visibility of the needle tip:

Dual Laser Marking  Embossed Marking 

HANDLES

2 types of handles available for individual comfort (only for single lumen flush & double lumen flush):

U shaped handle  Y shaped handle 



TRANSLUCENT TUBING

3 types of translucent tubing in any desired length (maximum length thick wall tubing is 120cm/47.2 inch):

Regular Tubing for standard oocyte aspiration set up
Soft Tubing allows control of the fluid
Thick Wall Tubing no kinking of the tubing

CORKS

2 types of silicone corks for a secure connection to test tubes:

Small	14x18x20 mm	0.55 x 0.71 x 0.79 inch		
Medium	21x27x30 mm	0.83 x 1.06 x 1.18 inch		



Contact our sales department for other or special enquiries!