

PROBET

ENDOMETRIAL HISTOLOGY

The Probet is an endometrial suction curette for a histologic lining biopsy.

FEATURES

- double headed piston secures a steady vacuum
- tissue collection through side hole
- smooth rounded distal tip and small size for minimal patient discomfort
- depth markings for ease in locating depth of insertion in the uterus



ORIGINAL PICTURE

GYNO INFO PROBET

Length	262 mm
Tube OD	2.93 mm
Sterile	ETO sterilized
Packaging	75/shipping box
Order number	#4164



GMP Belgium Probet 04/2011



Gynetics Medical Products n.v.
Rembergt Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine