

EMTRAC A

EMBRYO TRANSFER CATHETER WITH POSITIONER

The Emtrac A is a combination of a soft, flexible intrauterine catheter and a solid cervix catheter.


FEATURES

- positioner for depth adjustment
- a combination of a soft, flexible intrauterine catheter and a solid cervix catheter
- markings on outer tube for easy positioning
- handle with luer lock for a more secure fitting to syringes

ORIGINAL PICTURE

GYNO INFO EMTRAC A

Length	216 mm
Tube OD	1.60 mm
Total volume	0.42 ml
Total volume minus hub	0.258 ml
Sterile	sterilized by radiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	#4219 A



GMP Belgium Emtrac A 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine