

# EMTRAC DELPHIN

## EMBRYO TRANSFER CATHETER WITH POSITIONER AND SOFT UTERINE CATHETER

The combination of a soft flexible intrauterine catheter and a solid cervix catheter creates an universal catheter for easy insertion in the most difficult cervical tract.

### FEATURES

- positioner to prevent damage to fundus
- one piece catheter
- markings to indicate depth
- extra soft uterine catheter to ensure safe insertion through cervix
- handle with luer lock for a more secure fitting to syringes
- **available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance**

### TIP OPTIONS



CLASSIC

TIP-IN-SIGHT

### GYNO INFO EMTRAC DELPHIN

Length	210 mm
Tube OD	2.00 mm
Total volume	0.23 ml
Total volume minus hub	0.072 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

#### Order numbers

- Emtrac Delphin **#4219 Delphin**
- Emtrac Delphin Tip-in-Sight **#4219 Delphin TIS**

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.

