

GYNOSAMPLER

ENDOMETRIAL HISTOLOGY

The Gynosampler is the American name for the Probet which is sold in the USA.
The Gynosampler is an endometrial suction curette for a histologic lining biopsy.

FEATURES

- double headed piston secures a steady vacuum
- tissue collection through side hole
- smooth rounded distal tip and small size for minimal patient discomfort
- depth markings for ease in locating depth of insertion in the uterus



ORIGINAL PICTURE

GYNOSAMPLER

Length	262 mm
Tube OD	2.93 mm
Sterile	ETO sterilized
Packaging	250/shipping box (US only)
Order number	#9164 (US only)



GMP Belgium Gynosampler 04/2017