

PIPET-ID

MICROPIPETTOR

The Pipet-ID and Pipet-IDs are reusable IVF laboratory micropipettors developed to fit the "Strip-ID" micropipettes and most other popular micropipettes (Ø 0.9 mm OD)

FEATURES

- Available in two different sizes to fit every size of hand
 - **Pipet-ID** for regular to large hands
 - **Pipet-IDs** for small to regular hands
- Ergonomic design
- Adjustable aspiration volume (min ≈ 0.2 µl, max ≈ 3 µl)
- Total expulsion feature
- Stainless steel table holder included
- Spare parts kit included (includes extra seals, plunger wire and hexagon key)
- Designed for use with the "Strip-ID" micropipettes. Also fits most plastic micropipettes (Ø 0.9 mm OD, 0.6 mm ID)



ORIGINAL PICTURE

GYNO INFO PIPET-ID

Length (setting on position 3) 143 mm

Weight 57 gr

Order number #3000 Pipet-ID

GYNO INFO PIPET-IDs

Length (setting on position 3) 120 mm

Weight 43 gr

Order number #3000 Pipet-IDs



GMP Belgium Pipet-ID 10/2017



Colour differences may occur between batches due to variations in the colouring process and be considered as normal. GMP will not accept any claim or remark concerning these eventual differences.



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNÉTIQS[®]
Practise better Medicine