

# SEMTRAC SET W-240

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

### FEATURES

- soft bullet shaped tip to minimize trauma
- 5 **free** mandrels (B), to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter
- **available with Tip-in-Sight (TIS) technology for precise positioning of the catheter tip under ultrasound guidance**

### LOADING CATHETER TIP OPTIONS



### GYNO INFO SEMTRAC SET W-240

Guiding catheter (A) length	203 mm
tube OD	2.55 mm
Loading catheter (C) length	277 mm
tube OD	1.70 mm
total volume	0.26 ml
total volume minus hub	0.098 ml
Sterile	ETO sterilized
Packaging	75/shipping box
Quality control	MEA and LAL
Order numbers	
- Semtrac Set W-240	<b>#2000 W-240</b>
- Semtrac Set W-240 Trial	<b>#2000 W-240 T</b>
- Semtrac Set W-240 Tip-in-Sight	<b>#2000 W-240 TIS</b>

BEFORE USING REFER TO INSTRUCTIONS FOR USE (IFU) FOR INDICATIONS, CONTRA-INDICATIONS AS WELL AS WARNINGS AND PRECAUTIONS.

